




Control your protein level like magic!

ProteoTuner™ Protein Control System

기존 transcription 레벨에서의 제어(유전자 발현 유도 시스템)나 translation 레벨에서의 제어(RNAi)와는 완전히 다른 세포내의 표적 단백질의 양을 직접, 가역적으로 컨트롤 할 수 있는 획기적인 시스템!


- **Speed** 단백질의 ON/OFF는 물론, 발현량의 단계적인 제어가 가능
- **Reversibility** Ligand로 신속하고 가역적이며 정밀한 단백질 발현 조절
- **Simplicity** Vector 1종과 ligand (Shield1 또는 Guard1) 1종으로 구성된 간편한 시스템



* Ligand는 DD tag에 따라 Shield1과 Guard1 중 선택한다.

용도 표적 단백질의 기능 해석

원리 ProteoTuner™ System은 destabilization domain (DD)과 Shield1 또는 Guard1 중 1개의 ligand로 구성되어 있으며, 세포 내 목적 단백질의 양을 빠른 속도로 가역적인 조절이 가능한 시스템이다. 목적 유전자를 DD-tag나 DDG-tag가 있는 vector에 클로닝하여 발현시키면, tag된 단백질은 세포내 프로테오솜에 의해 분해되어 발현이 OFF 상태가 된다. 그러나 배지에 ligand (Shield1 또는 Guard1)를 넣으면 destabilization domain (DD)에 결합하여 목적 단백질이 분해되는 것을 막아 발현이 ON 상태가 된다. ProteoTuner™ Shield Systems과 ProteoTuner™ Guard Systems의 domain은 각기 다른 기술로부터 유래된 것으로 두 개의 시스템은 독립적으로 작용한다. 만약 두 시스템 모두 이용하면 2개의 단백질을 한꺼번에 조절할 수 있다.



Shield1

Shield1은 변형된 FKBP destabilization domain (DD)에 결합하여 DD-tag된 단백질을 프로테오솜 분해로부터 보호한다. 또한 Shield1 투입량의 조절을 통해 목적 단백질의 발현량을 조절 할 수 있다.

Guard1

Guard1은 변형된 *E. coli* DHFR destabilization domain (DDG)에 결합하여 DDG-tag된 단백질이 프로테오솜에 의해 분해되지 않도록 보호한다. 목적 단백질의 발현량을 Guard1으로 조절 할 수 있다.

※ 선택 가이드

ProteoTuner™ Guard는 transmembrane proteins에 최적화되어 있으며 일반적인 적용에는 ProteoTuner™ Shield를 권장한다.

Both Ligands	Shield1 Alone	Guard1 Alone	Neither Ligands

ProteoTuner™ System을 이용한 최신 논문

- 01 Agop-Nersesian, C. *et al.* (2010) *PLoS Pathog.* **6**(7): e1001029. Biogenesis of the inner membrane complex is dependent on vesicular transport by the alveolate specific GTPase Rab11B.
- 02 Agop-Nersesian, C. *et al.* (2009) *PLoS Pathog.* **5**(1): e1000270. Rab11A-controlled assembly of the inner membrane complex is required for completion of apicomplexan cytokinesis.
- 03 Almyog, G. and Nolan, G. P. (2009) *BioTechniques* **46**(1):44–50. Conditional protein stabilization via the small molecules Shld-1 and rapamycin increases the signal-to-noise ratio with Tet-inducible gene expression.
- 04 Armstrong, C. M. and Goldberg, D. E. (2007) *Nature Meth.* **4**(12):1007–1009. An FKBP destabilization domain modulates protein levels in *Plasmodium falciparum*.
- 05 Banaszynski, L. A. *et al.* (2006) *Cell* **126**(5):995–1004. A rapid, reversible, and tunable method to regulate protein function in living cells using synthetic small molecules.
- 06 Banaszynski, L. A. and Wandless, T. J. (2006) *Chem. Biol.* **13**(1):11–21. Conditional control of protein function.
- 07 Banaszynski, L. A. *et al.* (2008) *Nature Med.* **14**(10):1123–1127. Chemical control of protein stability and function in living mice.
- 08 Berdeaux, N. *et al.* (2007) *Nature Med.* **13**(5): 597–603. SIK1 is a class II HDAC kinase that promotes survival of skeletal myocytes.
- 09 Bhanot, H. *et al.* (2010) *Mol. Cancer Res.* [Epub ahead of print]. Induction of nonapoptotic cell death by activated Ras requires inverse regulation of Rac1 and Arf6.
- 10 Boulter, E. *et al.* (2010) *Nature Cell Biology.* **12**(5):477–483. Regulation of Rho GTPase crosstalk, degradation and activity by RhoGDI1.
- 11 Campeau, E. *et al.* (2009) *PLoS One* **4**(8):e6529. A versatile viral system for expression and depletion of proteins in mammalian cells.
- 12 Chu, B. W. *et al.* (2008) *Bioorg. Med. Chem. Lett.* **18**(22):5941–5944. Recent progress with FKBPderived destabilizing domains.
- 13 Daher, W. *et al.* (2010) *PLoS Pathog.* **6**(10): e1001132. Concerted action of two formins in gliding motility and host cell invasion by *Toxoplasma gondii*.
- 14 Dishinger, J. F. *et al.* (2010) *Nature Cell Biol.* **12**(7):703–710. Ciliary entry of the kinesin-2 motor KIF17 is regulated by importin-beta2 and RanGTP.
- 15 Doehn, U. *et al.* (2009) *Mol Cell.* **35**(4):511–522. RSK is a principal effector of the RAS-ERK pathway for eliciting a coordinate promotile/invasive gene program and phenotype in epithelial cells
- 16 Dvorin, J. D. *et al.* (2010) *Science.* **328**(5980):910–912. A plant-like kinase in *Plasmodium falciparum* regulates parasite egress from erythrocytes.
- 17 Evanko, Daniel (2006) *Nature Meth.* **3**(11):871. A protein alternative to RNAi.
- 18 Glass, M. *et al.* (2009) *Nat Methods.* **6**(8):577–579. *Conditional and reversible disruption of essential herpesvirus proteins.*
- 19 Grimley, J. S. *et al.* (2008) *Bioorg. Med. Chem. Lett.* **18**(2):759–761. Synthesis and analysis of stabilizing ligands for FKBP-derived destabilizing domains.
- 20 Haugwitz, M. *et al.* (2008) *BioTechniques* **44**(3):432–433. ProteoTuner: a novel system with rapid kinetics enables reversible control of protein levels in cells and organisms.

ProteoTuner™ Shield Systems

Product	Size	Cat. No.
ProteoTuner Shield System N	each	632172
ProteoTuner Shield System C	each	631072
ProteoTuner Shield System N (w/ AcGFP1)	each	632168
Lenti-X ProteoTuner Shield System N	each	632173
Lenti-X ProteoTuner Shield System C	each	631074
Lenti-X ProteoTuner Shield System N (w/ ZsGreen1)	each	632175
Retro-X ProteoTuner Shield System N	each	632171
Retro-X ProteoTuner Shield System N (w/ ZsGreen1)	each	632167
Shield1	60 µl	631037
	200 µl	631038
Shield1 (<i>in vivo</i>)	5 mg	632188
DD Monoclonal Antibody	50 µl	631073

ProteoTuner™ Guard Systems

Product	Size	Cat. No.
Lenti-X ProteoTuner Guard System N	each	631092
Lenti-X ProteoTuner Guard System C	each	631094
Retro-X ProteoTuner Guard System N	each	631093
Retro-X ProteoTuner Guard System C	each	631095
Guard1	60 µl	635051
	500 µl	635052

Takara

다카라코리아바이오테크놀로지(주)

고객지원센터 (학술/제품문의) 02 - 2081 - 2510
www.takara.co.kr support@takara.co.kr

본 사 대전지사
Tel. 02-2081-2510 Tel. 042-828-6525
Fax. 02-2081-2500 Fax. 042-828-6526