
SMARTer® Human TCR a/b Profiling Kit

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

SMART-Seq Technology

SMART = Switching Mechanism At 5' end of RNA Transcript

RACE Technology

RACE = Rapid Amplification of cDNA Ends

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


그림 1. Sequencing reads on target and correlation of clonotype count data for varying sample input amounts. A. Percentages of sequencing reads that map to CDR3 regions in either TCRα (blue) or TCRβ (purple) or that represent off target reads (gray). The protocol was performed on peripheral blood RNA: 10 ng, 100 ng, and 1,000 ng. For each RNA sample amount analyzed, ≥70% of sequencing reads mapped to a CDR3 region in either TCRα or TCRβ, with the 10 ng sample amount yielding the highest percentage of on target reads. B. Correlation of clonotype count data for varying sample input amounts. Another form of analysis involved plotting counts of the 100 most highly represented clonotypes for varying sample input amounts. Comparison of clonotype count data for the 100 ng and 1,000 ng sample amounts yielded a Pearson correlation coefficient of 0.80 and a Spearman coefficient of 0.80.

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그림 2. Assessing the reproducibility and sensitivity of the SMARTer approach.

A. The protocol was performed in replicate on PBMC RNA samples spiked at varying concentrations (10%, 1%, 0.1%, 0.01%, and 0.001%) with RNA obtained from a homogeneous population of Leukemic Jurkat T cells (TRBV12-3-TRBJ1-2 clonotype). The number of TRBV12-3-TRBJ1-2 specific sequence reads obtained for each spiked sample was normalized by subtracting the number of corresponding reads obtained for negative control samples consisting of unspiked PBMC RNA. Following Log_{10} transformation, plotting of the data (spike-in dilution vs. normalized read count) and linear regression analysis revealed a statistically significant correlation ($p = 3.93 \times 10^{-10}$, $R^2 = 0.99$) between the amount of spiked-in Jurkat RNA and the number of TRBV12-3-TRBJ1-2 specific sequence reads.

B. Count data, signal to noise ratios, and statistical analysis for TRBV12-3-TRBJ1-2 specific sequence reads obtained from spiked RNA samples. This result demonstrates that differences in the relative abundance of transcripts for a particular TCR clonotype are faithfully and reproducibly represented in sequencing libraries generated. Comparison of the number of TRBV12-3-TRBJ1-2 specific sequence reads obtained for the control vs. spike-in samples suggests that added Jurkat RNA at a concentration of 0.1% is detectable above background in the sequencing output ($p < 0.005$) at a depth of ~ 275,000 reads, evidence of the sensitivity afforded by the SMARTer approach.

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