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## What is microRNA?

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
**EXIQON**  
Seek Find Verify

## Exiqon microRNA 분석서비스

### Exiqon miRCURLY LNA™ microRNA Array, 7th gen

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**GENOCHECK**  
Genome Partner For Your Life

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
**그림 2. The structure of LNA.** The ribose ring is connected by a methylene bridge (orange) between the 2'-O and 4'-C atoms thus “locking” the ribose ring in the ideal conformation for Watson-Crick binding. When incorporated into a DNA or RNA oligonucleotide, LNA makes the pairing with the complementary strand more rapid and increases the stability of the resulting duplex.

	Same length, higher T <sub>m</sub>	Shorter length, similar T <sub>m</sub>
DNA	23-mer   tcgatcgattagctactagctga   T <sub>m</sub> : 60°C	23-mer   tcgatcgattagctactagctga   T <sub>m</sub> : 60°C
DNA/LNA™	23-mer   tcgatcgattAgctactagctga   T <sub>m</sub> : 64°C	16-mer   ---atcgattAgctActga---   T <sub>m</sub> : 49°C
DNA/LNA™	23-mer   tcgatcGattAgctAcGtaCgta   T <sub>m</sub> : 78°C	8-mer   .....aGctacGt.....   T <sub>m</sub> : 61°C


+ LNA™

DNA: atcg LNA: ATCG

**그림 3. Replace DNA with LNA for higher T<sub>m</sub>.** On the left, progressive substitution of DNA nucleotides with LNA increases the melting temperature of the oligonucleotide while maintaining the recognition sequence and specificity of the probe. On the right, LNA substitutions allow shortening of the probe while maintaining the same T<sub>m</sub>.

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**그림 4. The power of T<sub>m</sub> normalization.** The signal from DNA-based capture probes varies with GC content and results in poor detection of many microRNAs, whereas LNA probes offer robust detection of all microRNAs. Signal intensity from microarray experiments using LNA-enhanced (blue) or DNA-based (gray) capture probes. MicroRNA targets with varying GC content were added at 100amol each.

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
**그림 5. LNA microRNA inhibitors have high uniform potency.** The affinity of traditional full length microRNA inhibitors is highly influenced by the GC-content resulting in a Tm span of more than 40°C. In contrast, Exiqon's inhibitors span just 10°C around an optimal temperature.

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### DNA

- Real-time /quantitative PCR
- SNP detection/allele specific PCR
- Methylation analysis
- Bead-based applications
- Chromosomal FISH
- Comparative genome hybridization
- Proteomics of isolated chromatin segments (PICH)
- Antigene inhibition
- Mutagenesis

### mRNA

- Real-time /quantitative PCR
- Microarray analysis
- In situ hybridization
- Northern blotting
- Bead-based applications
- Fluorescence activated cell sorting
- Isolation
- Inhibition of RNA function
- RNA modification (frame shifting/exon skipping)
- DNAzymes

### ncRNA

- Real-time /quantitative PCR
- Microarray analysis
- In situ hybridization
- Northern blotting
- Fluorescence activated cell sorting
- Inhibition of RNA function
- RNA modification (frame shifting/exon skipping)

### miRNA

- Real-time /quantitative PCR
- Microarray analysis
- In situ hybridization
- Northern blotting
- Bead-based applications
- Inhibition of RNA function

#### PCR based approaches

#### Hybridization based approaches

#### In vivo based approaches

**그림 6. Proven LNA applications.** LNA is a powerful tool in many applications where standard DNA or RNA oligonucleotides do not have sufficient affinity or specificity. The figure shows an overview of some of the LNA applications that have been used for the study of RNA and DNA.

See how LNA™ works...  
Watch the LNA™ movie at  
[www.exiqon.com/e-talk](http://www.exiqon.com/e-talk)