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1. MicroRNA Isolation
2. Expression Analysis Array
3. Expression Analysis qPCR
4. Localization
5. Functional Analysis

1. MicroRNA Isolation (page 7)

miRCURY RNA Isolation Kits

- Get total RNA from a wide range of sources
- Ideal for microRNA research

2. Expression Analysis Array (page 8)

miRCURY LNA microRNA Array, 7th gen – hsa, mmu & mo

- Our sensitive and specific human, mouse and rat microRNA microarray

miRCURY LNA microRNA Hi-Power Labeling Kits

- Fast and simple labeling of total RNA
- Ideal for use with Exiqon's microRNA arrays

Spike-in microRNA Kits

- Improve the quality of your array data with these synthetic microRNAs

3. Expression Analysis qPCR (page 10)

Universal cDNA Synthesis & ExiLENT SYBR Green Master Mix kits

- Optimized reagents for use with the miRCURY LNA Universal RT microRNA PCR system

miRCURY LNA Universal RT microRNA PCR Primer Sets

- Individual microRNA PCR primer sets for quantification of microRNAs
- Design custom LNA primers for any small RNA online

Ready-to-Use PCR panels

- miRNome PCR panels
- microRNA Pick-&-Mix PCR panels
- Focus microRNA PCR panels

4. Localization (page 12)

miRCURY LNA microRNA Detection Probes

- Sensitive and specific probes for all microRNAs. Available with a wide variety of modifications

miRCURY LNA microRNA ISH Optimization Kits (FFPE)

- Kits for optimizing microRNA ISH from many sample sources
- First optimize your experiment, then use an LNA probe for your microRNA of interest

5. Functional Analysis (page 14)

miRCURY LNA microRNA Inhibitors

- Pre-designed and custom inhibitors for specific suppression of all microRNAs in miRBase and more

miRCURY LNA microRNA Power Inhibitors

- Our most potent inhibitors. Synthesized with PS backbones

miRCURY LNA microRNA Family Inhibitors

- Inhibitors of all microRNAs within a family

miRCURY LNA microRNA Target Site Blockers

- Unmatched high efficacy *in vitro* and *in vivo*
- Unrivalled performance - LNA TSBs do not catalyze RNase H-dependent mRNA degradation

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Sample Type	Recommended kit	RNA fraction	Recommended detection method
Serum Plasma Urine CSF Other biofluids	miRCURY™ RNA Isolation Kit - Biofluids	Small RNA (<1000 bp)	miRCURY LNA™ Universal RT microRNA PCR
Cultured cells Plant Tissues Brain and Adipose Tissue* Full Blood (non-coagulating)	miRCURY™ RNA Isolation Kit - Cell & Plant	Total RNA	miRCURY LNA™ Universal RT microRNA PCR or miRCURY LNA™ microRNA Arrays
Tissue	miRCURY™ RNA Isolation Kit - Tissue		

* For brain and adipose tissue extra Lysis Additive is required.

그림 1. Isolation kit selection guide. Select the right sample isolation kit and detection method for your samples.

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
그림 2. The miRCURY RNA Isolation Kit -Biofluids allows detection of more microRNAs compared to three other column based RNA isolation kits. RNA isolations were performed from the same plasma sample, using 4 different sample preparation kits. microRNA profiling was performed using the Serum/Plasma Focus microRNA PCR panel (168 microRNA assays in total). Call rate (percent of microRNAs detected) for each kit is shown. Only assays detected as present in all 5 replicates were calculated as a call.

표 1. Recommended starting material for different organisms /sample types

	Serum/plasma	Urine	CSF	Other biofluids*
Human samples	200 µl	200 µl	200 µl	200 µl
Rodent samples	50 µl**	200 µl	Not tested	50 µl**
RNA eluate for PCR***	4 µl	4 µl	8 µl	4 µl

*) For subsequent qPCR analysis, start volume may need to be adjusted to keep PCR inhibitors at a minimum.

**) Add RNase-free water to keep final volume at 200 µL

***) Volume used as input in a 20 µL RT reaction using Exiqon's Universal cDNA Synthesis Kit.

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
그림 3. High quality total RNA from difficult tissues. Total RNA was isolated from five mouse tissues using the miRCURY RNA Isolation Kit - Tissue. Brain and muscle tissues are usually considered difficult sources of RNA as they are very rich in lipids and fiber, respectively. However, great results were produced as seen from the average OD ratios (260/280: 1.85, 260/230: 1.7) and RIN values (8.8) of the five tissues.

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Organism name	Organism code	Common name	In miRBase 19	Array 7th gen
Homo sapiens	hsa	Human	2042	94%
Mus musculus	mmu	Mouse	1281	90%
Rattus norvegicus	mo	Rat	723	95%

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

그림 4. LNA-enhanced capture probes ensure robust detection of microRNAs. With DNA-based capture probes, half of microRNAs were either undetected or poorly detected. Signal strength (\log_2 signal/100amol target) from 660 synthetic microRNAs hybridized to Exiqon's microarray and Supplier A's DNA-based array are compared.


그림 5. T_m -normalized LNA capture probes. DNA capture probes (gray bars) for human microRNAs, have a T_m range of more than 30 °C and an average T_m of 64 °C (StDev 5 °C). The LNA capture probes (yellow bars) have a T_m range of only 10 °C and an average T_m of 71.5 °C (StDev=1.62 °C).

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
그림 6. The most sensitive array available. Due to optimally designed T_m normalized capture probes and extremely efficient labeling, the Exiqon array detects a significantly higher percentage of microRNAs than competitor arrays.

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
그림 7. Reliable microRNA expression profiles with as little as 30ng total RNA. Four different microarray experiments with varying amount of input RNA from oesophagus cancer (T) and normal adjacent (N) tissue were compared. A very high correlation is obtained when plotting the results from the experiment using 1000ng input RNA against those using 300, 100, and 30ng.

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표 2. Superior discrimination between microRNA family members. There is very little cross-hybridization between let-7 family members. The experiments were performed with synthetic let-7 spike-in microRNA (300 amol) in a background of tRNA.

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
그림 8. Great correlation between microarray and qPCR results. The array data was normalized (quantile normalization) and microRNAs with \log_2 ratios > 0.5 were included in the study. The qPCR data was normalized to reference genes. Only microRNAs that were detected (Cp < 36 for all replicates) were included. A total of 26 microRNAs were included in the study.

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
Step 1: First-strand synthesis (RT)

Step 2: Real-time PCR amplification


그림 9. Schematic outline of the miRCURY LNA Universal RT microRNA PCR System. A polyA tail is added to the mature microRNA template (step 1A). cDNA is synthesized using a PolyT primer with a 3' degenerate anchor and a 5' universal tag (step 1B). The cDNA template is then amplified using microRNA-specific and LNA-enhanced forward and reverse primers (step 2A). SYBR Green is used for detection (step 2B).

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
그림 11. Overview of the miRCURY LNA Universal RT PCR workflow. The PCR primer sets have been designed for optimal performance when used with Exiqon's ExiLent SYBR Green master mix. Use of other master mixes may affect the quality of the results. Ready-to-use (miRNome, Focus and Pick-&-Mix) panels can be replaced by individual PCR primers in this workflow.

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
그림 12. Accurate quantitation from 1pg total RNA starting material. Data from the amplification of 6 microRNAs in serial dilutions of human AM6000 total reference RNA are shown. All microRNA assays exhibit linear read-out with correlation coefficients $R(2) > 0.99$.

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
그림 13. Expression profiling of 742 microRNAs using 40 ng total RNA from tumor and normal FFPE sections. Real-time PCR was performed using triplicate RT reactions per sample on human miRNome microRNA panels I and II. Data from 424 microRNAs with Cp values <37 is included. Normalized expression is shown as fold changes in tumor compared to normal. Out of 424 microRNAs expressed, 144 were > 2-fold down-regulated in the tumor (blue dots) and 26 were > 2-fold up-regulated in the tumor (red dots).


그림 14. microRNA profiling in blood serum and plasma. Serum/Plasma Focus microRNA PCR Panels were used to profile 175 microRNAs commonly found in serum/plasma. MicroRNAs for sample quality control in orange and interesting microRNAs in red.

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표 3. Excellent discrimination between closely related microRNA family members. Examples of single nucleotide discrimination in the miR-181 family.

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
그림 15. Biomarker discovery workflow. Exiqon's qPCR system is designed for biomarker discovery, from screening in miRNome panels to final validation in a subset of samples by either Pick-&-Mix panels or individual assays.

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4. Localization

microRNA와 small RNA 검출을 위한 LNA-enhanced Detection system

miRCURY LNA microRNA Detection Probes for *in situ* hybridization & Northern Blotting

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
그림 16. LNA probes are superior to DNA probes. *A. thaliana* total RNA was hybridized with 32P-labeled DNA and LNA probes for miR-171. From Valoczi et al. 2004, Nucleic Acids Res. e175; reprinted with permission from Oxford University Press.

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그림 17. MicroRNA detection in zebrafish. Detection of miR-122a (top), miR-206 (middle) and miR-124a (bottom) using LNA probes in whole mount zebrafish embryos. Image kindly provided by Dr. Ronald Plasterk, Hubrecht Laboratory, The Netherlands.

그림 18. MicroRNA detection in chick. Specific detection of miR-206 in a Gallus gallus embryo using an LNA probe. miR-206 is detected in myotomal muscle cells (Ason et al. 2006).

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
그림 19. Double DIG labeling is more sensitive than single DIG labeling. hsa-miR-21 detection in tissue sections using an LNA probe with a double DIG (5' and 3') label at 40nM (A) or a single 3' DIG label at 80nM (B).

FFPE 샘플의 *in situ* hybridization 반응조건 최적화를 위한 miRCURY LNA microRNA ISH Optimization Kit(FFPE)

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
그림 20. Overview of the procedure. First, the tissue is "opened" using Proteinase K. In the hybridization step, the double DIG-labeled LNA probe binds specifically to its target microRNA. Alkaline phosphatase (AP)-conjugated anti-DIG antibodies are then added. This step is followed by NBT-BCIP development and optional counter-staining with Nuclear Red.

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
그림 21. miR-126 detection in colon wall. Kit 5 can be used to detect microRNAs in inflamed colon FFPE tissue. Here, it was used to detect miR-126. Staining was performed with NBT-BCIP (blue). Sections were counterstained with nuclear red.


그림 22. miR-145 detection in human colon. Kit 7 can be used to detect microRNAs in colon FFPE tissue. Here, miR-145 is detected in a human colon wall with underlying muscle layers. Staining was performed with NBT-BCIP (blue). Sections were counterstained with nuclear red.


그림 23. miR-205 detection in human breast carcinoma. Kit 8 can be used for detection of microRNAs in breast cancer FFPE tissue. Here, it was used to detect miR-205. Staining was performed with NBT-BCIP (blue). Sections were counterstained with nuclear red.

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표 4. Choosing the appropriate miRCURY LNA microRNA ISH Optimization Kit. The table indicates the tissue(s) in which each of the kits has been validated.

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miRCURY LNA microRNA Inhibitors & Power Inhibitors

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- Human library: 980 microRNA
- Mouse library: 739 microRNA
- Combined human and mouse library: 980 Human + 739 mouse microRNA

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
그림 24. Exiqon's microRNA inhibitors offer higher efficacy than leading competitors. MCF7 cells were cotransfected with a Firefly Luciferase-Reporter plasmid with a miR-21 target sequence, a transfection efficiency control and different concentrations of miR-21 inhibitors. Reporter gene expression was measured 48h after transfection with a Dual Luciferase assay.

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
그림 25. Silencing of common microRNAs. HeLa cells were transfected with a Firefly Luciferase-Reporter plasmid with a miR target sequence, a plasmid expressing Renilla luciferase (transfection efficiency control) and the corresponding miRCURY LNA microRNA Inhibitor in different concentrations. Reporter gene expression was measured 48 h after transfection with a Dual Luciferase assay. Ratios of Firefly and Renilla luciferase activity were calculated and normalized to values obtained with a Firefly Luciferase reporter with no miR target sequence (p-miR no target).


그림 26. Efficacy of microRNA Power Inhibitors. Cells transfected with miRCURY LNA microRNA Power Inhibitor against miR-181a (right panel) fail to differentiate normally (left panel). LHCN-M2 cells were plated and transfected with 50 nM of microRNA inhibitor or control. After 24 hours differentiation was induced by shifting to low serum medium. Seven days later the cells were fixed in formalin and permeabilized with Triton. Cells were stained with late differentiation marker (Myosin Heavy Chain (MHC)-Alexa488) and nuclei with Hoechst 33258. No MHC (green) is observed in miR-181a transfected cells. In addition, pictures without stains show that they also fail to form myotubes. Mir-181a is required to down-regulate HoxA11 which in turn is a repressor of MyoD. The absence of mir-181 causes a repression of MyoD, resulting in a failure of the cells to differentiate.

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
그림 27. LNA enhanced Target Site Blockers compete effectively with RISC for microRNA binding site. Target Site Blockers are antisense oligonucleotides designed to compete with microRNA/RISC by hybridizing to the microRNA target site of a particular mRNA.


그림 28. Target Site Blockers stimulate translation of specific mRNAs by masking the microRNA target sites. LNA enhanced high affinity target site blockers compete effectively with microRNA/RISC for the microRNA target site. In addition LNA distribution throughout the LNA/DNA mixmer ensures that the antisense oligonucleotide does not catalyze RNaseH dependent degradation of the mRNA. As a result the TSB will cause increased expression of the protein encoded by the targeted mRNA.

[License Notice www.Exiqon.com 참조]