LabelIT[®] Nucleic Acid Labeling Technology

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
그림 1. The LabelIT[®] Labeling covalently attaches labels to nucleic acids of all lengths. The LabelIT[®] chemistry covalently attaches fluorophores to the N7 position of guanines. The N7 position is not involved in base pairing, and the addition of the fluorophore to this position does not affect hybridization of labeled nucleic acids.

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
그림 2. Labeling density and fluorescent intensity are easily controlled using LabelIT[®] Reagents. Plasmid DNA was labeled with LabelIT[®] CX-Rhodamine Reagent at the indicated ratios (w:w) of labeling reagent to DNA and purified by ethanol precipitation. The density of labeling was determined by using two criteria: fluorescence intensity of the attached rhodamine and visible absorbance at 585 nm. A quenching effect is observed at the higher 5:1 ratio.

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■ **LabelIT[®] Tracker™ Intracellular Nucleic Acid Localization Kits**

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
그림 4. Tracking of plasmid localization and expression. COS-7 cells were transfected with LabelIT[®] Tracker™ Cy™5 labeled pEYFP-nuc and TransIT[®]-LT1 Transfection reagent in complete media. Images were acquired at 3, 8 and 20 hours post-transfection. The blue staining indicates the cellular localization of the labeled plasmid while the yellow signal is the expression of the nuclear yellow fluorescent protein (EYFP) reporter.

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