
Living Color® DsRed-Monomer Fluorescent Protein

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그림 1. Living Colors DsRed-Monomer is a monomeric protein. Panel A. Recombinant DsRed-Express and DsRed-Monomer (100 μ g) were analyzed by FPLC gel filtration chromatography. Overall absorption (A280) and chromophore excitation (A557) of the eluted material were monitored simultaneously. DsRed-Monomer elutes from the column at a retention time (39 min) corresponding to a molecular weight of 28 kDa. The calculated molecular weight of DsRed-Monomer is 26.8 kDa. DsRed-Express is a tetrameric protein that elutes at an earlier retention time (33 min) corresponding to a molecular weight of 89 kDa. Panel B. Pseudonative gel analysis of proteins. The oligomeric structures of proteins are preserved during SDS PAGE analysis if samples are kept at 4 °C and not boiled prior to loading on a gel. Boiled and unboiled recombinant proteins (7.5 μ g) were separated by SDS PAGE electrophoresis (12% acrylamide). In both the boiled (denatured) and unboiled (undenatured) samples, DsRed-Monomer and EGFP run as uniform bands of ~30 kDa due to their monomeric structures. The unboiled (undenatured) DsRed-Express runs at a much higher molecular weight than its boiled (denatured) counterpart due to its tetrameric structure.

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Protein	Excitation Max(nm)	Emission Max(nm)	Time to Detection (hr)	Relative Fluorescent Intensity	Quarternary Structure	Utility as a Reporter	Utility in Fusions
DsRed Monomer	556	586	12	Bright Fair	Monomer	Excellent	
DsRed Express	557	579	8-12	Extremely Bright	Tetramer	Excellent	Fair
DsRed2	563	582	24	Extremely Bright	Tetramer	Good	Fair

표 1. Spectral properties of DsRed proteins

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표 2. Filter sets for DsRed-Monomer detection


그림 3. Cells expressing DsRed-Monomer are easily detected by standard BD FACS™ analysis. HEK 293 cells were transfected with pDsRed-Monomer-N1 and analyzed using the 568 nm laser line of the BD FACSVantage™ SE Flow Cytometry System 24 hr post-transfection.

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
그림 4. DsRed-Monomer is soluble when expressed in mammalian cells, HeLa cells were transfected with pDsRed-Monomer-N1 and fixed in 4% paraformaldehyde 24 hr post-transfection, DsRed-Monomer displays an even, consistent, and homogeneous distribution.


그림 5. DsRed-Monomer performs well in protein fusions and multicolor applications. Various DsRed-Monomer fusions were transiently transfected in HeLa cells and visualized by fluorescence microscopy. Cells were fixed in 4% paraformaldehyde 24 hr post-transfection. Panel A, DsRed-Monomer-Caveolin fusion, Panel B, DsRed-Monomer-Rab5 GDP/GTP GTPase exchange factor homologue fusion, Panel C, DsRed-Monomer-Golgi (trans Golgi stack) and AcGFP1-Nuc (nucleus), Panel D, DsRed-Monomer-Lamin B fusion and AcGFP1-Actin

Reference

1. Reef Coral Fluorescent Protein Vectors(July 2003) Clontechiques XVIII(3):6-7.
2. Matz, M. V., et al. (1999) Nature Biotechnol, 17:969-973, Erratum in: Nature Biotechnol, (1999) 17:1227.
3. Living Colors™ DsRed-Express(July 2003) Clontechiques XVII(3):16.
4. Living Colors™ DsRed Polyclonal Antibody (Jan 2003) Clontechiques XVIII(1):11.

■ Product list

Living Colors® DsRed-Monomer	Products Size	Cat. No.
pDsRed-Monomer Vector	20 µg	632467
pDsRed-Monomer-N1 Vector	20 µg	632465
pDsRed-Monomer-C1 Vector	20 µg	632466
pDsRed-Monomer-Actin Vector	20 µg	632479
pDsRed-Monomer-Golgi Vector	20 µg	632480
pDsRed-Monomer-F Vector	20 µg	632493
pDsRed-Monomer-Hyg-N1 Vector	20 µg	632494
pDsRed-Monomer-Hyg-C1 Vector	20 µg	632495
pDsRed-Monomer-N In-Fusion Ready Vector	1 µg	632498
pDsRed-Monomer-C In-Fusion Ready Vector	1 µg	632499
pRetroQ-DsRed Monomer-N1 Vector	20 µg	632507
pRetroQ-DsRed Monomer-C1 Vector	20 µg	632508
pDsRed-Monomer-Mem Vector	20 µg	632512
pDsRed-Monomer-Mem Hyg Vector	20 µg	632513
pDsRed-Monomer-F Hyg Vector	20 µg	632514
rDsRed-Monomer Protein	100 µg	632503