PreMix WST-1 Cell Proliferation Assay System

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TaKaRa Code MK400 2,500회 (25 ml x 1 bottle of PreMix WST-1)

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
그림 2. Determination of the activity of Human interleukin-2 (IL-2) on the mouse T cell line CTLL-2. Mouse CTLL-2T Cells were seeded in 96 well microtiter plates (tissue culture grade, 96 wells, flat bottom) at a concentration of 4x10³ cells/well, with 100 ul of culture medium containing various amounts of IL-2 (final concentration e.g. 0.005-10 ng/ml). Incubate cells for 48 hrs at 37℃ and 5%, and then 10 ul/well PreMix WST-1 was added. After an additional incubation for 4 hrs under the same condition, A₄₅₀ was measured in a multiwell plate reader.


그림 3. Measurement of cytotoxic activity of human TNF- α to WEHI-164 cells. Mouse fibrosarcoma cell line WEHI-164 cells were preincubated at a concentration of 1×10^6 cells/ml in culture medium with actinomycin C1 for 3 hrs at 37°C and 5% CO₂. Cells were seeded at a concentration of 5×10^4 cells/well in 100 μ l culture medium containing actinomycin C1 and various amounts of TNF- α into microtiter plate. After incubation for 24 hrs at 37°C and 5% CO₂, 10 μ l/well PreMix WST-1 was added and then cells were incubated an additional 4 hrs under the same conditions. The absorbance was measured at A₄₅₀ in a multiwell plate reader.

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1. LDH Cytotoxicity Detection Kit MK401 2,000회

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