Ready-To-Glow™ Secreted Luciferase Reporter Systems

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
그림 1. High signal intensity and stability using secreted *Metridia* luciferase. CHO cells were plated into 96-well plates and transiently transfected with CMV-driven constructs encoding non-secreted firefly luciferase, non-secreted *Renilla* luciferase, and sequence-optimized secreted *Metridia* luciferase, 24hr after transfection, luciferase activity in equivalent sample amounts was analyzed by addition of the recommended substrate. The signal was measured at different timepoints over a period of 45 min. (neg = negative control.)



그림2. Ready-To-Glow™ Secreted Luciferase Reporter System workflow comparison with other non-secreted luciferase reporter systems.

Ready-To-Glow™ Dual Secreted Reporter System

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Comparable Reporter Activity is suitable for time course assays

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

그림3. Similar secretion kinetics of *Metridia* secreted luciferase and SEAP enable accurate comparisons of the relative timing of promoter activity.

HeLa cells were plated into 6-well plates and transiently transfected with either the pMetLuc-Control vector or the pSEAP-Control vector. Media samples from the transfected cells were collected from 3 wells at each time point by removing enough media to run either the luciferase or the SEAP assay. Each sample was tested in triplicate, using a white bottom 96-well microtiter plate on a Turner BioSystems Veritas® Luminometer

Monitor two promoters simultaneously

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Retain live cells for other applications

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
그림 4. Monitoring activation of two promoters simultaneously. HEK 293 cells co-transfected with pNFkB-TA-MetLuc and pCRE-SEAP constructs were treated with fresh media alone or with media containing either 1,000 ng/ml TNF- α or 10 μ M forskolin. Culture supernatant samples were collected 7 hr later and assayed as in Figure 1, using a BD Monolight™ 3096 Luminometer.

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Product	TaKaRa Code	Size
Ready-To-Glow Secreted Luciferase pMetLuc Vector Kit	631729	each
Ready-To-Glow Secreted Luciferase Reporter Assay	631726	100 rxns
	631727	500 rxns
	631728	1,000 rxns
Ready-To-Glow Secreted Luciferase Reporter System	631730	100 rxns
	631731	500 rxns
	631732	1,000 rxns
Ready-To-Glow Dual Secreted Reporter Assay	631734	500 rxns
Ready-To-Glow Dual Secreted Reporter Vector Kit	631735	4 x 20 ug

TALON®

Bacterial Expression and Purification Kits

- ▶ Clone, express, and purify – all with a single Kit
- ▶ Based on the well-characterized pET system