In-Fusion™ 2.0 PCR Cloning Kits

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- Clones PCR fragments up to 12 kb in just minutes
- No PCR fragment purification or restriction digestion required
- Use with Advantage HD DNA Polymerase for superior accuracy and efficiency

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Advantage HD DNA Polymerase is used to generate an amplified PCR fragment containing your gene of interest (GOI). Following PCR amplification, a 30 min treatment with our proprietary Cloning Enhancer (37°C for 15 min, then 80°C for 15 min) eliminates the need to purify the PCR product using spin columns. After treatment, the heat-inactivated PCR product is ready for a 30 min In-Fusion cloning reaction (37°C for 15 min, then 50°C for 15 min). The lyophilized format requires only the addition of PCR product and linearized vector prior to incubation. During the In-Fusion cloning reaction, the In-Fusion enzyme produces single-stranded regions at the ends of the PCR fragment and vector that are fused together due to the 15 bases of homology. The resulting reaction can be used to transform *E. coli* without further treatment. This simple method is well-suited for cloning a GOI contained within a large PCR fragment.

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A 12 kb PCR fragment was amplified from Adeno-X™ using Advantage HD DNA Polymerase. Following a proprietary Cloning Enhancer treatment, this fragment was combined in an In-Fusion cloning reaction with pDNR-Dual (4.9 kb) that had been prelinearized by digestion with *Sal* and *HindIII*. The In-Fusion reaction was then transformed into Fusion-Blue™ Competent Cells, and transformants were selected on LB/Amp plates. The resulting colonies were picked at random and screened by colony PCR, by amplifying a 2 kb fragment within the expected 12 kb fragment. Eight out of eight, or 100%, of the randomly picked clones contained the desired recombinant plasmid. Lane M : 1 kb DNA ladder molecular weight marker; Lanes 1~8: sampled colonies.

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
그림 3. In-Fusion 2.0 PCR cloning eliminates the need for spin column purification and provides superior cloning efficiency.

A 2 kb DNA fragment was amplified from Adeno-X using Advantage HD DNA Polymerase, and cloned into pDNR-Dual using the following three methods: Panel A. The amplified 2 kb fragment was spin column-purified and then cloned into pDNR-Dual via an In-Fusion Dry-Down reaction. Panel B. The amplified 2 kb fragment was purified using our proprietary Cloning Enhancer and then cloned into pDNR-Dual via an In-Fusion Dry-Down reaction. Panel C. The amplified 2 kb fragment was cloned into pDNR-Dual via an In-Fusion Dry-Down reaction without any prior treatment. Equivalent volumes of each In-Fusion reaction were transformed into Fusion-Blue Competent Cells and transformants were selected on LB/Amp plates. The Cloning Enhancer-treated sample (Panel B) yields significantly more recombinant clones (fivefold more) than that obtained using spin column purification (Panel A), which shows improvement over the untreated sample (Panel C).

제품명	Size	TaKaRa Code
In-Fusion 2.0 Dry-Down PCR Cloning Kit	8 rxns	639609
In-Fusion 2.0 CF Dry-Down PCR Cloning Kit*	24 rxns	639607
	96 rxns	639608
Fusion-Blue Competent Cells	24 rxns	636700
	96 rxns	636758
Advantage HD DNA Polymerase	200 rxns	639241

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- pDNR-CMV (linearized control)
- 1.1 kb Unpurified Control Insert
- Cloning Enhancer
- In-Fusion Dry-Down Reactions (8 rxns, 24 rxns, or 96-well plate)
- Fusion-Blue™ Competent Cells (8 rxn size only)