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- Strong and tightly regulated suppression of any target gene
- Save weeks of time by creating a stable, inducible cell line in only one step
- Ideal for cases where gene suppression may be lethal

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
그림 1. The Knockout Single Vector Inducible RNAi System uses a modified form of the tightly regulated, tetracycline-controlled gene expression system. In the absence of Dox, tTS binds the tetO sequences within the TRE/U6 promoter and actively silences transcription of the shRNA. When Dox is added to the culture medium, tTS dissociates from the TRE, relieving transcriptional suppression and allowing high level transcription of the shRNA from the hybrid TRE/U6 promoter.

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■ Effective & Sensitive Gene Suppression

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
그림 2. Doxycycline-induced knockdown of lamin A/C in HeLa cells. Panel A, A stable HeLa cell line that expresses an anti-lamin A/C shRNA was produced using the Knockout Single Vector System. Cells were grown in the absence or presence of Dox (1 µg/ml) for the times indicated. They were then harvested and analyzed by Western blot using β-actin expression as a control. Suppression of lamin A/C expression is evident after 6 hr of treatment, and is virtually abolished after 72 hr. Lane 1: control, Lanes 2-4: 6 hr, 48 hr, and 72 hr, respectively. Panel B, The anti-lamin A/C HeLa cell line was grown in the absence or presence of increasing concentrations of Dox for 72 hr, then harvested and analyzed by Western blot. Lane 1: parental cells. Lane 2-4: 0, 0.1, and 1 µg/ml Dox, respectively.


그림 3. Sensitive doxycycline-induced knockdown of luciferase activity.  
HEK 293 cells were transiently co-transfected with the pSingle-tTS-shRNA vector expressing an anti-luciferase shRNA and a pCMV-luciferase expression vector at a ratio of 1:1. Cells were grown in the absence or presence of 0–1 ug/ml Dox for 72 hr, then harvested and lysed to measure luciferase activity. Luciferase activity was reduced by 67% at 1 ng/ml Dox and by 88% at 1 ug/ml.

#### Reference

1. Liu, C.Y., et al. (1996) *Genes Devel.* 10(14):1835–1843.
2. Gossen, M., & Bujard, H. (1992) *Proc. Natl. Acad. Sci. USA* 89(12):5547–5551.
3. Freundlieb, S., et al. (1999) *J. Gene Med.* 1(1):4–12.
4. Witzgall, R., et al. (1994) *Proc. Natl. Acad. Sci. USA* 91(10):4514–4518.
5. Wiznerowicz, M., & Trono, D. (1999) *J. Virol.* 77(16):8957–8961.
6. pTRE-Tight Vectors (April 2003) *Clontechiques XVIII*(3):13–14.
7. Yao, F., et al. (1998) *Hum. Gene Ther.* 9(13):1939–1950.

**Tet-On® & Tet-Off®  
Advanced**

- ▶ Inducible RNAi and cDNA expression
- ▶ Up to 10,000-fold inducibility
- ▶ Highly sensitive and controllable