효율적인 Transfection Method

Calcium Phosphate Transfection of Neuronal Cells

Min Jiang & Gong Chen
 Department of Biology, Huck Institutes of Life Sciences,
 The Pennsylvania State University, University Park, Pennsylvania 16802.
 Correspondence should be addressed to G.C. (gongchen@psu.edu).

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Forming a Homogeneous DNA-Ca²⁺ - Phosphate Precipitate

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
그림 2. Formation and subsequent dissolution of the DNA-Ca²⁺-phosphate precipitate. Panel A, Continuous vortexing when mixing DNA with Ca²⁺ and phosphate buffer results in large clusters of precipitate (examined after a 1 hr incubation). Panels B & C, Formation of an optimal DNA-Ca²⁺-phosphate precipitate through gentle vortexing during mixing (image taken after 1 hr incubation). Panel D, Dissolution of the precipitate with slightly acidic transfection medium pre-equilibrated in a 10% CO₂ incubator (Reproduced from references 1 and 12).

Dissolving the DNA-Ca²⁺ Phosphate Precipitate

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
그림 3. High transfection efficiency is achieved in low-density hippocampal cultures with our improved protocol. Panels A-C, Phase-contrast micrographs. Panels D-F, Fluorescent images of GFP-transfected cells in three independent transfections. Note that the majority of neurons in the local field (microislands) are transfected. Neurons were cultured for 10-15 days according to previously described methods^(11, 12) Scale bar, 50 μm (Reproduced from references 1 and 12).

High Transfection Efficiency in Low-Density Hippocampal Cultures

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Transfection of Mature Neurons in Culture

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
그림 4. Successful transfection of both mature and immature neurons. Panels A & B. Transfection of young neurons with EGFP at 2 and 7 days in vitro (d.i.v.). Panels C & D. Transfection of mature neurons with EGFP at 24 and 82 d.i.v. Scale bar, 50 mm (Reproduced from reference 1).

High Transfection Efficiency in Neurons but not Glial Cells

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
그림 5. Transfection efficiency in glial cells is very low. Panel A. Only a few glial cells (arrowhead) are transfected in neuronal cultures with a monolayer of astrocytes. Panel B. In Banker-type cultures where neurons usually do not contact glial cells directly, transfection efficiency remains high (21.6%). Panels C & D. Many constructs can be successfully transfected using our protocol. Illustrated here are EGFP-Rab3a (Panel C) and EGFP-dynamin (Panel D) (Reproduced from reference 1).

Product	Size	TaKaRa Code
CalPhos™ Mammalian Transfection Kit	each	631312

References

- Jiang, M. & Chen, G. (2006) *Nat. Protoc.* 1(2):695-700.
- Xia, Z. *et al.* (1996) *J. Neurosci.* 16(17):5425-5436.
- Micheva, K. D. *et al.* (2003) *Nat. Neurosci.* 6 (9):925-932.
- Chen, W. G. *et al.* (2003) *J. Neurosci.* 23 (7):2572-2581.
- Holz, R. W. *et al.* (2000) *J. Biol. Chem.* 275 (23):17878-17885.
- Passafaro, M. *et al.* (2003) *Nature* 424 (6949):677-681.
- Goetze, B. *et al.* (2004) *J. Neurobiol.* 60(4):517-525.
- Craig, A. M. (1998) *Transfecting Cultured Neurons*, MIT Press, Cambridge, MA.
- Washbourne, P. & McAllister, A. K. (2002) *Curr. Opin. Neurobiol.* 12(5):566-573.
- Chen, Y. *et al.* (2003) *Cell* 115(1):37-48.
- Chen, G. *et al.* (2004) *Proc. Natl. Acad. Sci. USA* 101(4):1063-1068.
- Jiang, M. *et al.* (2004) *Gene Ther.* 11(17):1303-1311.
- Kohrmann, M. *et al.* (1999) *J. Neurosci. Res.* 58(6): 831-835.
- Watanabe, S. Y. *et al.* (1999) *Neurosci. Res.* 33(1):71-78.