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
그림 1. The Great EscAPe™ SEAP Chemiluminescence Kit 2.0 is more sensitive and has greater signal stability than the original kit. HEK 293 cells were transiently transfected with a vector encoding SEAP under the control of the constitutive CMV promoter. The media supernatant was collected after 48 hr from transfected and mock-transfected cells. Samples were assayed with the original and new kits and protocols, and analyzed on a Turner BioSystems Modulus™ Microplate Multimode Reader.

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Live Cell Assay Provides Flexibility

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
그림 2. SEAP activity is highly specific with a > 20-fold increase in activity. HEK 293 cells were transiently transfected with a promoter construct containing CRE driving the expression of SEAP, or mock-transfected. 12 hr posttransfection, the cell culture media was replaced by media containing either 10 μ M forskolin or plain cell culture media and incubated for 7 hr. Forskolin causes an increase in the level of cytosolic cAMP, which in turn activates CRE, driving the expression of SEAP, which is detectable in the media culture supernatant. The samples were assayed using the new kit and protocol and analyzed on a BD Monolight™ 3096 Luminometer.

Product	Size	TaKaRa Code
Great EscAPe™ SEAP	50 rxns	631736
Chemiluminescence Kit 2.0	300 rxns	631737
	1,000 rxns	631738

Components

- 5 x Dilution Buffer
- SEAP Substrate Solution
- Positive Control Placental AP
- Protocol-at-a-Glance (PT3954-2)

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- Ready-To-Glow™ Dual Secreted Reporter Assay (TaKaRa Code 631734)
- Ready-To-Glow™ Dual Secreted Reporter Vector Kit (TaKaRa Code 631735)
- Fusion-Blue™ Competent Cells (TaKaRa Code 636700 & 636758)
- Supercharge EZ10 Competent Cells (TaKaRa Code 636756)

Reference

1. Ready-To-Glow™ Dual Secreted Reporter System (October 2006) *Clontechiques* XXI (3):1.