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
그림 1. Living Colors® DsRed-Monomer is a monomeric protein, Panel A, Recombinant DsRed-Express and DsRed-Monomer (100 μg) were analyzed by FPLC gel filtration chromatography. Overall absorption (A280) and chromophore excitation (A557) of the eluted material were monitored simultaneously. DsRed-Monomer elutes from the column at a retention time (39 min) corresponding to a molecular weight of 28 kDa. The calculated molecular weight of DsRed-Monomer is 26.8 kDa, DsRed-Express is a tetrameric protein that elutes at an earlier retention time (33 min) corresponding to a molecular weight of 89 kDa, Panel B, Pseudonative gel analysis of proteins. The oligomeric structures of proteins are preserved during SDS PAGE analysis if samples are kept at 4 °C and not boiled prior to loading on a gel. Boiled and unboiled recombinant proteins (7.5 μg) were separated by SDS PAGE electrophoresis (12% acrylamide). In both the boiled (denatured) and unboiled (undenatured) samples, DsRed-Monomer and EGFP run as uniform bands of ~30 kDa due to their monomeric structures. The unboiled (undenatured) DsRed-Express runs at a much higher molecular weight than its boiled (denatured) counterpart due to its tetrameric structure.

Protein	Excitation Max (nm)	Emission Max (nm)	Time to Detection (hr)	Relative Fluorescent Intensity	Quaternary Structure	Utility as a Reporter	Utility in Fusions
DsRed-Monomer	556	586	12	Bright	Monomer	Fair	Excellent
DsRed-Express	557	579	8-12	Extremely Bright	Tetramer	Excellent	Fair
DsRed2	563	582	24	Extremely Bright	Tetramer	Good	Fair

표 1. Spectral properties of DsRed proteins

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표 2. Filter sets for DsRed-Monomer detection

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
그림 3. Cells expressing DsRed-Monomer are easily detected by standard BD FACS™ analysis. HEK 293 cells were transfected with pDsRed-Monomer-N1 and analyzed using the 568 nm laser line of the BD FACSVantage™ SE Flow Cytometry System 24 hr post-transfection.

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그림 4. DsRed-Monomer is soluble when expressed in mammalian cells. HeLa cells were transfected with pDsRed-Monomer-N1 and fixed in 4% paraformaldehyde 24 hr post-transfection. DsRed-Monomer displays an even, consistent, and homogeneous distribution.

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
그림 5. DsRed-Monomer performs well in protein fusions and multicolor applications. Various DsRed-Monomer fusions were transiently transfected in HeLa cells and visualized by fluorescence microscopy. Cells were fixed in 4% paraformaldehyde 24 hr post-transfection. Panel A, DsRed-Monomer-Caveolin fusion, Panel B, DsRed-Monomer-Rab5 GDP/GTP GTPase exchange factor homologue fusion, Panel C, DsRed-Monomer-Golgi (trans Golgi stack) and AcGFP1-Nuc (nucleus), Panel D, DsRed-Monomer-Lamin B fusion and AcGFP1-Actin.

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