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■ TALON[®] Resin Offers 6xHis Protein Purification Under Native Conditions Using β -mercaptoethanol as a Reducing Agent

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High Protein Purity and No Metal Leakage with β -Mercaptoethanol

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Better Yields Using β -ME with TALON than Ni-NTA

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
그림 1. TALON 6xHis DHFR purification under native conditions with increasing concentrations of β -ME. N-terminal 6 x His-tagged mouse DHFR (19.5 kDa) was expressed in *E. coli*. Cell cultures (25 ml) were disrupted in 2 ml of sonication buffer (pH 8.0) by freezing and thawing. Clarified lysate (2.66 ml at 4 mg/ml total protein) was applied to 200 μ l of TALON in a gravity flow column, pre-equilibrated with the sonication buffer. The 6xHis protein/resin complexes were washed 3 x with sonication buffer, pH 8.0. All bound protein was eluted with 600 μ l of 100 mM EDTA. The indicated concentrations of β -ME were included in the buffer. Samples were electrophoresed on a 12% polyacrylamide/SDS gel, and the gel was stained with Coomassie blue. Lane 1: 20 μ l of cell lysate, Lanes 2, 4, 6, & 8: 20 μ l of flowthrough. Lanes 3, 5, 7, & 9: 5 μ l of eluant.


그림 2. Yields of recombinant 6xHis DHFR from cell extracts purified by TALON versus Nickel-NTA in the presence of β -ME. N-terminal 6xHis DHFR was expressed and purified under native conditions as described in 그림 1. Protein concentrations were determined by Bradford assay[®]. Yields are expressed as a percentage of total protein in the cell lysate.

References

1. TALON Resin (April 1995) *Clontechiques X*(2):8-9.
2. Porath, J. (1992) *Protein Express, Pur.* 3(4):263-281.
3. Hochuli, E., et al. (1987) *J. Chromat.* 411:177-184.
4. Sawadogo, M., et al. (1995) *Genet. Eng.* 17:53-65.
5. Bradford, M. (1976) *Anal. Biochem* 72:248.

Rapid One-Step IMAC Purification Using TALON[®] Superflow

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
그림 3. Rapid one-step purification of 6xHis-GFPuv on TALON Superflow. Whole cell extract containing 6xHis-GFPuv was loaded in 50 mM sodium phosphate, 0.3 M NaCl, pH 8.0 (arrow 1), washed in this buffer plus 10 mM imidazole (arrow 2), and eluted in the same buffer plus 150 mM imidazole (arrow 3).

References

1. Porath, J., et al. (1975) *Nature* 258(5536): 598-599.
2. Hochuli, E., et al. (1987) *J. Chrom.* 411: 177-184.
3. TALON Metal Affinity Resin (April 1995) *Clontechiques X* (2):8-9.
4. Kain, S. R. (1997) *Biotechnol. Intl.* 1:79-84.


그림 4. Analysis of TALON Superflowpurified fractions. SDS-PAGE analysis was performed on the fractions represented in 그림 1. Lane 1 : extract, Lane 2 : nonadsorbed material (peak I), Lane 3 : wash (peak II), Lane 4 : eluate (peak III), Lane M : MW standard.

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