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
그림 2. Overview of developing the stable cell lines of the Tet-Off Advanced System.

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
그림 3. Gene induction in the Tet-Off Advanced and Tet-On Advanced Systems. The Tet-controlled transactivators for these systems (tTA-Advanced and rTA-Advanced) are fusion proteins derived from a wild-type or mutant Tet repressor (TetR), respectively. Each DNA-binding TetR domain is joined to three minimal transcription activation domains of the herpes simplex virus VP16 protein, and each has been optimized for expression in mammalian cells. For Tet-Off Advanced, the uninduced basal state is maintained in the presence of doxycycline, which prevents the tTA-Advanced protein from binding to TREmod sequences in pTRE-Tight. Removal of doxycycline permits tight binding and induces high-level transcription. In exact contrast, the Tet-On Advanced System is activated by doxycycline.

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#### • Tet System Approved Fetal Bovine Serum

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#### 【Tet Advanced start-up system】

- Tet-Off® Advanced Inducible Gene Expression System 1 Kit 630934
- Tet-On® Advanced Inducible Gene Expression System 1 Kit 630930

#### 【premade Tet-On® Advanced cell line】

- HEK293 Tet-On® Advanced Cell Line 1 ml 630931
- HepG2 Tet-On® Advanced Cell Line 1 ml 630932
- MCF7 Tet-On® Advanced Cell Line 1 ml 632108
- HeLa Tet-On® Advanced Cell Line 1 ml 632110
- HeLa Tet-Off® Advanced Cell Line 1 ml 632111
- MCF7 Tet-Off® Advanced Cell Line 1 ml 632109
- HEK293 Tet-Off® Advanced Cell Line 1 ml 632107
- Hep G2 Tet-Off® Advanced Cell Line 1 ml 632106

#### 【Tet System Approved Fetal Bovine Serum】

- Tet System Approved FBS, USDA-Approved 50 ml 631107
- Tet System Approved FBS, US-Sourced 50 ml 631105

#### 【TetR 검출용 monoclonal antibody】

- TetR Monoclonal Antibody 40 µg 631108

## Inducible MicroRNA and Living Colors® Fluorescent Protein Reporters Using the Tet-On® Advanced System

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
그림 6. Organization of the anti-lamin A/C and anti-Luc (L1) miRNA expression cassettes. Anti-luciferase and anti-lamin A/C miRNA-like sequences, based on the miR-30 miRNA<sup>®</sup>, were joined to the coding regions of AcGFP1 or DsRed2 cDNAs in the pTRE-Tight expression plasmid. Sense (S) and anti-sense (AS) target sequences are connected by a 19 nt loop sequence. Transcripts are processed to yield either an active miRNA or a mature mRNA for the fluorescent protein.

### A Tet-On Advanced System for miRNA Delivery

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### Linking miRNAs to Fluorescent Reporters

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
그림 7. Doxycycline-induced suppression of lamin A/C in HeLa cell lines. Stable, inducible HeLa cell lines were produced that expressed the Tet-On Advanced transactivator and either the pTRE-Tight-lamin A/C-AcGFP1 or -DsRed2 construct. Cells were grown in the presence or absence of Dox and assayed after 72 hr by Western Blot analysis. Lane 1 : Parental cells. Lane 2 : Anti-lamin A/C-AcGFP1 cells without Dox. Lane 3 : Anti-lamin A/C-AcGFP1 with Dox (1 µg/ml). Lane 4 : Anti-lamin A/C-DsRed2 without Dox. Lane 5 : Anti-lamin A/C-DsRed2 with Dox (1 µg/ml).

### Inhibition of Endogenous HeLa Cell Lamin Expression

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
그림 8. Doxycycline-induced fluorescence in anti-lamin A/C HeLa cell lines. The HeLa cell lines described and treated as in 그림 7 were visualized by fluorescence microscopy.  
 Panel A, Anti-lamin A/C-AcGFP1 cells without Dox.  
 Panel B, Anti-lamin A/C-AcGFP1 cells with Dox (1 µg/ml).  
 Panel C, Anti-lamin A/C-DsRed2 cells without Dox.  
 Panel D, Anti-lamin A/C-DsRed2 cells with Dox (1 µg/ml).  
 Dox treatment-induced strong fluorescence in cells that also exhibited lamin A/C knockdown.

### Suppression of a Luciferase Transgene

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
그림 9. Doxycycline treatment reduces luciferase activity in HEK-293 cells transfected with the Anti-Luc-DsRed2 reporter.  
 HEK-293 cell lines stably transformed with Tet-On Advanced were transiently cotransfected with pCMV-Luc and pTRE-Tight-L1-DsRed2 in triplicate at a ratio of 1:2, and then grown in the presence of Dox (1 µg/ml) for 72 hr. The L1 miRNA induced by Dox treatment caused a 85% reduction in luciferase activity. S.D. = standard deviation.


그림 10. Visualization of HEK-293 cells expressing the Anti-Luc-DsRed2 reporter in the absence and presence of doxycycline.  
 Panel A, In the absence of Dox, the Anti-Luc-DsRed2 reporter is undetectable by fluorescence microscopy.  
 Panel B, In the presence of Dox (1 µg/ml), Anti-Luc-DsRed2 expression reaches high levels.

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