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In-Fusion™ 2.0 Cloning System

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
그림 1. Comparison of Clontech, Invitrogen, and Stratagene PCR cloning methods. Panel A, Clontech's In-Fusion™ 2.0 cloning method may be used with any linearized vector. The PCR product is amplified with primers that contain 15 bp extensions homologous to the ends of the linearized vector. Panel B, pcDNA3.1 Directional TOPO Expression Kit cloning method from Invitrogen. The PCR product is amplified with a forward primer containing a 5' CACC extension. Panel C, Stratagene's StrataClone Ultra™ Blunt PCR Cloning Kit method. No primer extensions are required and cloning is not directional. The linear vector arms, in combination with the insert, are circularized by Cre recombinase after transformation into the StrataClone Competent Cells. The Invitrogen pcDNA3.1 and Stratagene vector arms are charged with topoisomerase I.

Benefits of the In-Fusion™ Cloning System

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Superior Cloning Efficiency for a Range of Insert Lengths

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
그림 2. PCR fragments amplified for comparison cloning. Fragments of 0.5 kb, 1 kb, 2 kb, 4 kb, and 6 kb were amplified using Advantage® HD polymerase. Panel A, PCR inserts generated for the In-Fusion™ 2.0 and Stratagene cloning. (Note that these inserts include the 15 bp homology extensions for In-Fusion cloning into pDNR-Dual. These extensions have no impact on Stratagene's cloning method.) Panel B, PCR inserts generated for Invitrogen cloning with the included CACC addition to the forward primer for directional cloning. 25 ng of linearized pAdeno-X was used as template for these PCR products. 1 µl of a 10 µM solution of each primer was used per 50 µl PCR reaction. Lanes M and M' : 1 kb ladder molecular weight markers; Lanes 1 and 1' : 0.5 kb fragment; Lanes 2 and 2' : 1 kb fragment; Lanes 3 and 3' : 2 kb fragment; Lanes 4 and 4' : 4 kb fragment; Lanes 5 and 5' : 6 kb fragment. The StrataClone Ultra method has no particular primer design criteria and thus any PCR product could be cloned. For this method, we chose to clone the same PCR product that was generated for the In-Fusion Cloning method.

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표 1 : Transformation Efficiency (Number of Colonies)¹

PCR Insert Size (kb)	In-Fusion™ 2.0 (1/10 Transformation Plated)	Stratagene (1/10 Transformation Plated)	Invitrogen (Whole Transformation Plated)
0,5	~ 950	19	236
1,0	~1,600	11	164
2,0	~1,570	20	~540
4,0	~1,050	18	298
6,0	~1,480	57	101

¹ The numbers of colonies generated for the indicated inserts are shown. Note that for the In-Fusion™ 2.0 and the StrataClone products, only 1/10 of the total transformation mixes were plated. For the Invitrogen product, the entire transformation was plated.

표 2 : Ratio of Plasmids Containing the Correct Insert¹

PCR Insert Size (kb)	In-Fusion™ 2.0		Stratagene		Invitrogen	
	Ratio	%	Ratio	%	Ratio	%
0,5	8/8	100	4/8	50	6/8	75
1,0	8/8	100	0/8	0	4/8	50
2,0	8/8	100	4/8	50	0/8	0
4,0	8/8	100	0/8	0	6/8	75
6,0	4/8	50	3/8	37,5	0/8	0

¹ The ratio and percent of positive inserts relative to the number of clones sampled for each of the three products are shown. The orientation of the clones for the nondirectional StrataClone method was not determined.


그림 3. Identification of clones containing the correct insert. The results for the 4 kb insert are shown here as representative. Eight randomly selected clones from each cloning method were digested with the restriction enzymes *Hind* III and *Sal* I to determine if the correct insert was present, (*Hind* III and *Sal* I sites were introduced in the 5' and 3' PCR primers.)

Panel A, Clontech's In-Fusion™ 2.0 Cloning System,

Panel B, Stratagene's StrataClone cloning system,

Panel C, Invitrogen's Directional Topo cloning system,

The ratios of positive clones were as follows : In-Fusion™ 2.0 : 8/8, StrataClone : 0/8, Invitrogen : 6/8, Lane M : 1 kb ladder molecular weight markers.

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In-Fusion™ 2.0 Dry- Down PCR Cloning kit Components

- pDNR-CMV (linearized control)
- 1,1kb Unpurified Control Insert
- Cloning Enhancer
- In-Fusion™ Dry-Down Reactions (8 rxns, 24 rxns, or 96-well plate)
- Fusion-Blue™ Competent Cells (8 rxn size only)

Product	Size	Cat.No
In-Fusion™ 2.0 Dry-Down PCR Cloning Kit	8 rxns	639609
In-Fusion™ 2.0 CF Dry-Down PCR Cloning Kit	24 rxns	639607
	96 rxns	639608
Advantage® HD Polymerase Mix	200 rxns	639241

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