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Cell Biology Group, Clontech Laboratories, Inc., in collaboration with Molecular Devices, now a part of MDS Analytical Technologies

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NF-κB Translocation은 Proteasome-의존적

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
그림 1, NF-κB activation by TNF-α requires proteasomal activity. High levels of *Metridia* luciferase signal were only observed in the absence of ALLN and the presence of TNF-α. When the proteasome is inactivated by ALLN (monitored by increasing levels of ZsProSensor-1 fluorescence), the NF-κB signaling pathway cannot respond effectively to TNF-α stimulation. ALLN = N-acetyl-leucyl-leucyl-norleucinal.

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
그림 2, Requirement for active proteasomes in TNF-α-induced, NF-κB-dependent signaling. Inactive NF-κB is sequestered in the cytoplasm by I-κB. I-κB must be phosphorylated upon TNF-α induction (1) and degraded by the proteasome (2) in order for NF-κB to translocate to the nucleus and initiate signaling (3). Alternatively, when the proteasome is inhibited by the peptide ALLN (4), I-κB is not degraded and NF-κB cannot translocate (5). The status of the proteasome (active or inactive) can be monitored based on ZsProSensor-1 levels.

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1. Adams, J. (2002) *Trends Mol. Med.* **8** (4 Suppl):S49-S54.
2. Haas, M, et al. (1998) *J. Leukocyte Biol.* **63** (3):395-404.
3. Sun, L, and Carpenter, G. (1998) *Oncogene*, **16** (16):2095-2102.
4. Proteasome Sensor Vector. (April 2003) *Clontechiques XVII* (2):14.
5. Matz, M, V, et al. (1999) *Nature Biotech* **17** (10):969-973. *Erratum in Nature Biotech*, **17** (12):1227.
6. Li, X, et al. (1998) *J. Biol. Chem.* **273** (52):34970-34975.
7. Living Colors® Cell Lines, (April 2004) *Clontechiques XVII* (2):14.
8. Measurement of proteasome inhibition in live cells in Molecular Devices microplate fluorometers, (2006) *SpectraMax Application Note*. Molecular Devices.
9. Ready-To-Glow™ Secreted Luciferase System, (July 2006) *Clontechiques XXI* (2):12-13.
10. Ready-To-Glow™ Dual Secreted Reporter System, (October 2006) *Clontechiques XXI* (3):1.
11. Detection Limit and Linear Range of Ready-To-Glow Secreted *Metridia* Luciferase Determined on the Molecular Devices SpectraMax® L, (July 2007) *Clontechiques XXII* (3):20-21.
12. SpectraMax® M5/M5e multi-detection reader: A five-mode microplate reader with three-mode cuvette port, *Molecular Devices Data Sheet* (2006).
13. SoftMax® Pro software : The industry standard in microplate data analysis, *Molecular Devices Data Sheet* (2006).

Ready-To-Glow™ Secreted Reporter Systems - FAQs

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1. Detection Limit and Linear Range of Ready-To-Glow™ Secreted *Metridia* Luciferase Determined on the Molecular Devices SpectraMax L. (July 2007) *Clontechiques* **XXI** (3):20-21.
2. Ready-To-Glow™ Secreted Luciferase System.(July 2006) *Clontechiques* **XXI** (2):12-13.
3. Ready-To-Glow™ Dual Secreted Reporter System.(October 2006) *Clontechiques* **XXI** (3):1.

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제품명	Size	TaKaRa Code
Proteasome Sensor Vector	20 µg	632425
Living Colors HEK 293 ZsGreen Proteasome Sensor Cell Line	1 vial	631535
Ready-To-Glow Secreted Luciferase pMetLuc Vector Kit	each	631729
Ready-To-Glow Secreted Luciferase Reporter Assay	100 rxns	631726
	500 rxns	631727
	1,000 rxns	631728
Ready-To-Glow Secreted Luciferase Reporter System ¹	100 rxns	631730
	500 rxns	631731
	1,000 rxns	631732
Ready-To-Glow Dual Secreted Reporter Vector Kit	4 x 20 µg	631735
Ready-To-Glow Dual Secreted Reporter Assay	500 rxns	631734

¹The Reporter System combines the Reporter Assay and the pMetLuc Vector Kit components.

※ License Notice : [21, 26, 27]