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Cell Biology Group  
Clontech Laboratories, Inc.

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
그림 1. Shield1 does not affect viability, even at 5–10 times the recommended concentration. Each indicated cell type ( $1 \times 10^4$  cells/well of a 96-well plate) was treated with 0, 1, 3, or 5 μM Shield1. Staurosporine (50 nM), a known apoptosis inducer, was used as a negative control. 48 hr later, cell proliferation was measured using the Premixed WST-1 Cell Proliferation Reagent.


그림 2. Shield1 has a stabilizing effect on secreted *Metridia* luciferase protein, implying a similar mechanism for Shield1-based stabilization of secreted and cytosolic proteins. Shield1 stabilizes secreted proteins in a tunable manner. Cells were transiently transfected with the DD-*Metridia* luciferase construct, and treated with increasing concentrations of Shield1. Six hr later, media samples were collected and the amount of secreted DD-*Metridia* luciferase in the media was determined using the Ready-To-Glow Secreted Luciferase Reporter Assay. RLU = relative light units.

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
그림 3. Tunable rescue of the ZAP70 signal transduction pathway in the ZAP70 deficient 116 Jurkat cell line.


그림 4. Rescue of the ZAP70-dependent signaling pathway with the ProteoTuner system. The level of rescue is proportional to the amount of ZAP70 that is rescued, 116 Jurkat cells (ZAP70-negative) were infected or mock-infected with a construct for DD-ZAP70, cultured in varying concentrations of Shield1, and induced with PMA and CD3. The extent of recovery was measured by IL-2 secretion into the culture media, via an ELISA assay.

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\* The number of reactions depends on the concentration of Shield1 used. At the maximum

\* License Notice: [1], [2], [3], [5], [6], [7], [8], [9], [13], [14], [16]

# Berthold Detection Systems

## Your Partner for Luminescence Measurement

### FB12 Tube Luminometer

The Gold Standard for Glow Luminescence

- Ultrasensitive Photon Counting Detector
- Total Flexibility in Sample Format
- Standalone Operation or PC Operation

### Sirius Tube Luminometer

Highest Performance and Flexibility  
same features as FB12, PLUS:

- Built-in Printer
- Up to Two Automated Reagent Injectors
- Safe Operation through Sample Presence Detection

### Orion L Microplate Luminometer

All You Need for Luminescence  
Measurement

- Ultrasensitive Photon Counting Detector
- Up to Two Automated Reagent Injectors
- Sophisticated PC Software

### Orion II Microplate Luminometer

The Top Instrument in the Field  
same features as Orion L, PLUS:

- Up to Four Automatic Reagent Injectors
- Two Temperature Control Options
- And Many More Powerful Functions