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
그림 1. HUVEC cells treated with camptothecin for 24h. Cells were assayed using the ToxiLight assay kit according to the manufacturer's instructions. The data produced was compared to a 100% lysed control sample and is expressed as a % lysis. The data is representative of the mean of 3 separate experiments (containing six replicates per dose) ±SEM


그림 2. ATP levels in HUVEC cells treated with camptothecin for 24h. Cells were assayed using the ViaLight Plus bioassay kit according to the manufacturer's instructions. The data is representative of the mean of 3 separate experiments (containing six replicates per dose) ±SEM

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
그림 3. Comparison of ViaLight Plus and ToxiLight Assays using HUVECS dosed with camptothecin. The ATP levels indicated by the ViaLight Plus RLUS reduce steadily in a dose dependent manner. At the lower drug doses the AK released from the cells is relatively low compared with that of the control only increasing dramatically at the highest drug doses.


그림 4. Live dead staining of HUVEC cells dosed with Camptothecin. A, 0 nM camptothecin; B, 500 nM camptothecin; C, 1000 nM camptothecin; D, 5000 nM camptothecin. The top level of photographs shows the calcein staining. The number of cells fluorescing green decreases as the level of drug increases. Conversely in the lower photographs the red stain of cells with compromised membranes increases with increasing drug.

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제품명	용량	TaKaRa Code
ViaLigh Plus Assay Kit	500 회	LT07-221
	1,000 회	LT07-121
	10,000 회	LT07-321
	500 회 (TC plate 5개 포함)	LT17-221
ToxiLight Non-Destructive Cytotoxicity BioAssay Kit	1,000 회	LT07-117
	500 회	LT07-217
	500 회 (TC plate 5개 포함)	LT17-217
ToxiLight 100% Lysis Control Set	100 ml (200 tests)	LT07-517

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Lonza Nottingham, Ltd