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그림 1. The HisTALON Cartridges (1 ml) provide highly efficient and specific method for purifying his-tagged proteins.

표 1. Comparison of HisTALON Cartridge (1 ml) Protein Yield & Activity to that of Competitor Cartridge<sup>1</sup>

Vendor	Starting Sample		Flowthrough		Wash		Eluate	
	Protein (mg)	AcGFP1 (RFU)	Protein (mg)	AcGFP1 (RFU)	Protein (mg)	AcGFP1 (RFU)	Protein (mg)	AcGFP1 (RFU)
Clontech	22.3	202,167	21.1	9,363	0.1	292	1.1	198,114
Vendor G <sup>2</sup>	12.7	82,845	12.2	4,980	0.1	215	0.6	81,111
Vendor Q <sup>2</sup>	12.4	88,300	11.3	4,560	0.2	3,600	0.7	78,320
Vendor P <sup>3</sup>	20.6	201,174	18.6	10,304	0.14	653	1,311	193,214

<sup>1</sup> Extraction and chromatography were performed according to the respective vendor's recommendations(see the Figure 2 caption for the Clontech sample extraction procedure).

<sup>2</sup> Extraction in presence of lysozyme.

<sup>3</sup> Extraction with the product-specific recommended extraction buffer.

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
그림 2. HisTALON Cartridges provide highly reproducible and rapid his-tagged protein purification. Protein extraction and chromatography were performed according to the User Manual, with 0.4 g pellets of E. coli cells expressing his-tagged AcGFP1. These cells were extracted in the recommended volumes of extraction buffer and centrifuged at 10,000 x g for 20 min, at 4°C. For each run, extract from 250 mg of cells was consequently applied to the cartridge using the recommended loading buffer and flow rate, on an FPLC system from GE Healthcare (used for all the purifications described in this article). After the wash and elution steps were completed, all fractions, as well as the original samples, were analyzed for protein content<sup>(7)</sup> and relative fluorescence (from active AcGFP1) on a 96-well fluorescence spectrophotometer. A high level of reproducibility was observed between runs, with differences of less than 5% in both protein content and fluorescent signal (data not shown). The entire process, including extraction and purification, was completed in less than 1 hour.


그림 3. HisTALON Cartridges yield higher purity than a competitor. Fractions from the purification of his-tagged AcGFP1 on respective HisTALON and Vendor P cartridges were compared via electrophoresis on a 4-15% SDS polyacrylamide gel stained with Coomassie Blue. Lanes M: Bio-Rad Precision Plus Protein All Blue Standards. Lane 1: HisTALON extract, Lane 2: HisTALON flowthrough, Lane 3: HisTALON eluate. Lane 4: Vendor P extract, Lane 5: Vendor P flowthrough, Lane 6: Vendor P eluate. A contaminating high molecular weight band appeared in Lane 6 (Vendor P eluate), but not in Lane 3 (HisTALON eluate), indicating that the HisTALON eluate was of higher purity—consistent with the data in Table I.

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
그림 4. HisTALON Cartridges effectively purify mammalian samples. A 200 mg HEK 293 cell pellet sample was extracted with our TALON xTractor Buffer. The extract, spiked with 25 µg of AcGFP1 and containing more than 10 mg of mammalian proteins from the HEK 293 cells, was loaded on a HisTALON cartridge. Fractions were analyzed via electrophoresis on a 4-15% SDS polyacrylamide gel and visualized with a silver staining kit (both from Bio-Rad Laboratories). Lanes M: Bio-Rad Precision Plus Protein All Blue Standards, Lane 1: HEK 293 extract spiked with 25 µg of AcGFP1, Lane 2: flowthrough, Lane 3: eluate. The single band in Lane 3 demonstrates the high purity of the eluate.

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제품명	용량	TaKaRa Code
HisTALON Cartridge	5 x 1ml	635650
HisTALON Buffer Set	20 purifications	635651
HisTALON Cartridge Purification Kit <sup>1</sup>	5 x 1ml	635649

<sup>1</sup> The HisTALON Cartridge Purification Kit combines the HisTALON Cartridge (containing 5 x 1 ml cartridge) and the HisTALON Buffer Set.

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1. Porath, J. *et al.* (1975) *Nature* 258 (5536) : 598-599.
2. Hochuli, E. *et al.* (1987) *J. Chrom.* 411:177-184.
3. Hochuli, E. *et al.* (1988) *Bio/Technology* 6(11):1321-1325.
4. Chaga, G. *et al.* (1994) *Protein Eng.* 7(9):1115-1119.
5. Froelich, C. J. *et al.* (1996) *Biochem. Biophys. Res. Commun.* 229(1):44-49.
6. Stephens, L. R. *et al.* (1997) *Cell* 89 (1): 105-114.
7. Bradford, M. M. (1976) *Anal. Biochem.* 72:248-254.

\* License Notice: [18]