Yeast Two-Hybrid 시스템의 결정판 Matchmaker™ Gold Yeast Two-Hybrid System

[1] Matchmaker Gold Yeast Two-Hybrid System

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
그림 1. Yeast two-hybrid system design. Library-derived, transcription-activating prey fusion proteins that interact with the DNA-binding bait fusion protein activate the expression of reporter genes.


그림 2. Four reporters give Matchmaker Gold its high stringency. Interacting bait and prey fusion proteins drive the expression of four different reporters from three different GAL4-responsive promoters (M1, G1, and G2), which are stably integrated in the genome of the reporter strain, Y2HGold. Aureobasidin A (AbA) resistance and the two auxotrophic reporters for histidine and adenine biosynthesis confer growth selection in the presence of AbA and on histidine- and adenine-deficient media, while the α -galactosidase reporter produces blue colonies in the presence of X- α -Gal.

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continued...


그림 3. The Mate & Plate Protocol. To screen a Matchmaker Mate & Plate Library, an aliquot of the library in the Y187 strain (*MAT α*) is simply mixed with a bait-expressing culture of the Y2HGold strain (*MAT α*). The mated strains are cultured overnight and plated on selective agar medium (e.g. -Leu/-Trp + AbA + X-alpha-Gal).

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
그림 4. Secondary Matchmaker Gold screening confirms high numbers of positive clones. A Y2HGold bait containing the POU domain from the mouse Oct4 transcription factor (BD-POU_{Oct4}) was used to screen the Mate & Plate Universal Mouse (Normalized) Library for Oct4-binding proteins. 32 colonies from a low stringency primary screen (DDO + AbA, 60 ng/ml + X-alpha-Gal) were selected and replated/patched onto fresh low stringency medium (Panel A) and onto high stringency medium (QDO + AbA, 60 ng/ml + X-alpha-Gal) (Panel B) to confirm the expression of all four Matchmaker Gold reporters. Of the 32 originally selected colonies, 25 were confirmed positive for the 4 reporters. DDO = Double dropout medium: SD/-Leu/-Trp. QDO = Quadruple dropout medium: SD/-Ade/-His/-Leu/-Trp.

표 1. Mouse Oct4-Binding Proteins Identified in Matchmaker Gold Library Screening¹

Protein	Function
E21/Ubc9	Small ubiquitin-related modifier (SUMO) enzyme involved in sumoylation of Oct4; regulates Oct4 stability
PIAS1	An E3 ligase protein inhibitor of activated STAT1; a potent inhibitor of Oct4-mediated transcriptional activation
PSMB5	Proteasome beta5 subunit; may mediate the interaction of Oct4 with proteasomes, which regulate cellular processes through protein degradation.

¹ Mate & Plate Library – Universal Mouse (Normalized).

Yeast Two-Hybrid 시스템의 결정판 Matchmaker™ Gold Yeast Two-Hybrid System

continued...


그림 5. Titration of Aureobasidin A for three different Matchmaker protein pairs. Y2HGold yeast clones co-expressing one of three bait/prey fusion protein pairs: a negative control (BD-POU_{mOct4} + AD-null); a pair of putative interactors (BD-POU_{mOct4} + AD-E21); or a positive control (BD-p53 + AD-T-Ag), were grown on DDO agar (SD/-Trp/-Leu) in the presence of increasing concentrations of AbA. The relative ability of each strain to grow in the presence of AbA reflects the strength of the interactions between the bait and prey proteins. Libraries are usually screened in the presence of AbA at 60 ng/ml, which demonstrates definitive selection for interacting hybrids. POU_{mOct4} = POU domain from the mouse transcription factor, Oct4.

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제품명	용량	TaKaRa Code
Matchmaker Gold Yeast Two-Hybrid System	each	630489
Make Your Own "Mate & Plate" Library System	5 rxns	63049

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Matchmaker Gold Yeast Two-Hybrid System

- pGBKT7 DNA-BD Cloning Vector
- pGADT7 AD Cloning Vector
- pGBKT7-53 Control Vector
- pGBKT7-Lam Control Vector
- pGADT7-T Control Vector
- Y2HGold Yeast Strain
- Y187 Yeast Strain
- YPDA Broth
- YPDA with Agar
- SD/-Trp with Agar
- SD/-Leu with Agar
- Yeastmaker Yeast Transformation System 2

continued...

[2] Mate & Plate Yeast Two-Hybrid cDNA Libraries

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
그림 6. The Mate & Plate Protocol. To screen a Mate & Plate Library, an aliquot of the library in the MAT α Y187 strain is simply mixed with a bait-expressing, MAT α reporter strain culture (Y2HGold or AH109). The two strains are co-cultured overnight and then plated on selective agar medium.


그림 7. Normalization reduces the abundance of cDNAs derived from highly expressed genes. Universal cDNA, synthesized using RNA from mixed human tissues, was analyzed before and after normalization on a NimbleGen Homo sapiens microarray (Cat. No. A4542-00-01). Data are shown for the 7,000+ genes which exhibited greater than 10,000 Mean Fluorescence Units (MFU). The signal intensities for approximately 3,300 of the most highly represented genes were significantly reduced following normalization, reflecting a preferential reduction of these abundant cDNAs.

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
그림 8. DSN-Normalization removes highly abundant cDNAs. Normalized (Lanes N) and non-normalized (Lanes C) human HeLa S3 cDNAs were compared using virtual Northern blot analysis and 32P-labeled probes. The levels of these two highly abundant cDNAs were sharply reduced following normalization.

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제품명	용량	TaKaRa Code
Mate & Plate Library - Universal Human (Normalized)	2 x 1 ml	630481
	5 x 1 ml	630480
Mate & Plate Library - Universal Mouse (Normalized)	2 x 1 ml	630482
	5 x 1 ml	630483
Mate & Plate Library - HeLa S3 (Normalized)	5 x 1 ml	630479
Mate & Plate Library - Human Bone Marrow	5 x 1 ml	630477
Mate & Plate Library - Human Fetal Brain	5 x 1 ml	630469
Mate & Plate Library - Human Heart	5 x 1 ml	630471
Mate & Plate Library - Human Liver	5 x 1 ml	630468
Mate & Plate Library - Human Skeletal Muscle	5 x 1 ml	630473
Mate & Plate Library - Human Testis	5 x 1 ml	630470
Mate & Plate Library - Human Ovary	5 x 1 ml	630474
Mate & Plate Library - Mouse Embryo 11-day	5 x 1 ml	630478
Mate & Plate Library - Mouse Embryo 17-day	5 x 1 ml	630476

continued...

[3] Convenient Premixed Yeast Two-Hybrid Media Sets

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그림 9. Specialized yeast media mixes are available in convenient, ready-to-go pouches.

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Yeast Media Set 2 (0.5 L pouches)

- YPDA Broth (2 each)
- YPDA with Agar
- SD/-Leu Broth
- SD/-Leu with Agar
- SD/-Trp Broth
- SD/-Trp with Agar
- SD/-Leu/-Trp (DDO; 10 each)
- SD/-Ade/-His/-Leu/-Trp (QDO)

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Yeast Media Set 2 Plus

- Yeast Media Set 2
- Aureobasidin A (1 mg)
- X-alpha-Gal (250 mg)

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1. Takesako, K. *et al.* (1991) *J. Antibiot.* 44 (9):919-924.
2. Fields, S & Song, O. (1989) *Nature* 340 (6230):245-246.
3. Chien, C. T. *et al.* (1991) *Proc. Nat. Acad. Sci. USA* 88(21):9578-9582
4. Pretransformed Normalized Matchmaker™ Library (January 2007) *Clontechiques* XXII (1):21-23.
5. Pretransformed Normalized Matchmaker™ Libraries (January 2008) *Clontechiques* XXIII (1):14-16.
6. High-Performance Reference RNA and cDNA (July 2008) *Clontechiques* XXIII (2):18-20.
7. Matchmaker™ Gold Yeast Two-Hybrid System. (January 2009) *Clontechiques* XXIV (1):1-3.

제품명	용량	TaKaRa Code
Yeast Media Set 2	each	630494
Yeast Media Set 2 Plus	Each	630495
YPDA Broth	10 pouches	630306
YPDA with Agar	10 pouches	630307
SD/-Leu Broth	10 pouches	630310
SD/-Leu with Agar	10 pouches	630311
SD/-Trp Broth	10 pouches	630308
SD/-Trp with Agar	10 pouches	630309
SD/-Leu/-Trp Broth	10 pouches	630316
SD/-Leu/-Trp with Agar	10 pouches	630317
SD/-Ade/-His/-Leu/-Trp Broth	10 pouches	630322
SD/-Ade/-His/-Leu/-Trp with Agar	10 pouches	630323
Aureobasidin A	1 mg	630466
X-alpha-Gal	100 mg	630462
	250 mg	630463

* License Notice: [12], [15], [17]