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Cell Biology Group
Clontech Laboratories, Inc.

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
그림 1. Dynamic polymerization and depolymerization of actin. Self-assembly of actin filaments occurs at the plus end of an existing actin filament, as monomeric actin is incorporated. Conversely, disassembly occurs at the minus end where actin monomers depolymerize from the filament, causing a continuous rearrangement of the actin filament network.

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
그림 2. Experimental design. In the absence of Shield1 only mCherry-Actin is present to be incorporated into newly forming actin filaments. When Shield1 is added, DD-AcGFP1-Actin is stabilized, and is therefore incorporated together with mCherry-Actin into newly formed actin filaments.

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continued...


그림 3. Controlling actin stability with the ProteoTuner system. In the absence of Shield1, DD-AcGFP1-Actin is not present (Panels B & F) despite a normal, mCherry-labeled actin filament network (Panels A & E). In the presence of Shield1, DD-AcGFP1-Actin is stabilized, and present in the actin network (Panels D & H) in addition to mCherry-labeled actin (Panels C & G). Cells were fixed using 4% paraformaldehyde and imaged with a 40X objective on a Zeiss Axioskop microscope using the green Chroma filter set HQ460/40, Q490LP, HQ515/30 and the red Chroma filter set HQ 540/40X, 570DCPL, D600/50M.

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* The number of reactions depends on the concentration of Shield1 used. At the maximum suggested concentration (1,000 nM) 60 µl = 30-plus reactions and 200 µl = 1,000-plus reactions in a six-well plate.

* Shield1 (60 µl) is included with all ProteoTuner Systems.

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