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
그림 1. SMARTer cDNA synthesis occurs in a single-step reverse transcription reaction. Following amplification, SMARTer cDNA can be used for a variety of downstream applications.

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
그림 2. Typical gel profile of ds cDNA synthesized with the SMARTer PCR cDNA Synthesis Kit (Panel A) and the SMARTer Pico PCR cDNA Synthesis Kit (Panel B), using the Control Human Placental Total RNA as a template. 2 ng (Panel A) or 1 ng (Panel B) of the Control Human Placental Total RNA was subjected to first-strand cDNA synthesis and purification, followed by PCR amplification as described in the user manual for each kit. 5 μl of each PCR product was electrophoresed on a 1.2% agarose/EtBr gel in 1X TAE buffer following the indicated number of PCR cycles. Lanes M : 1 kb DNA ladder size markers. The arrows indicate the strong band at 900 bp typically seen for human placental total RNA.

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제품명	용량	TaKaRa Code
SMARTer PCR cDNA Synthesis Kit	10 rxns	634925
	20 rxns	634926
SMARTer Pico PCR cDNA Synthesis Kit	10 rxns	634928
SMARTer RACE cDNA Amplification Kit	10 rxns	634923
	20 rxns	634924
Universal Primer Mix	100 rxns	634922

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
그림 3. The SMARTer II A Oligo increases the overall yield of the 5'-RACE PCR products. Note the difference in intensity between the 2.6 kb 5'-RACE product in Lane 2 (produced by the SMART II A oligo) and that in Lane 6 (produced by the SMARTer II A oligo); both bands resulted from the amplification of only 2 ng of human placental total RNA template. In the assay shown, first-strand cDNA synthesis and 5'-RACE PCR amplification were performed as described in the user manual (PT4096-1). Lanes 1 & 5: No RNA 5'-RACE control, Lanes 2 & 6: Transferrin Receptor (TFR) 5'-RACE starting from 2 ng total RNA. Lanes 3 & 7: TFR 5'-RACE starting from 50 ng total RNA. Lanes 4 & 8: 5'-cDNA internal control using 2 ng total RNA as template. Lane M: 1 kb DNA Ladder.

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- Chenik, A. et al. (1998) In *Gene Cloning and Analysis by RT-PCR*. Eds Siebert, P. & Larrick, J. (BioTechniques Books, MA) Ch 22.
- Be SMART™ About First-Strand cDNA Synthesis (January 2009) *Clontechiques* XXIV(1): 15 - 17.

* License Notice : [16]