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
그림 2. Cloning options when using the EP and CEP Starter Kits. Features of the multiple cloning sites of our pEcoli Expression Vectors are shown (not to scale). Convenient restriction sites (indicated) allow inclusion or removal of tags as needed. The unshaded region is removed in our pEcoli Linear Expression Vectors using *Sal* I and *Hind* III restriction enzymes.

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
그림 3. Excellent purity is obtained when using HisTALON Cartridges for purification. The CEP Starter Kit was used to generate BL21 *E. coli* expressing N-terminal 6xHN-tagged AcGFP1, as seen in Panels A & B; and C-terminal 6xHN-tagged AcGFP1, as seen in Panels C & D. HisTALON Cartridges were used to purify the expressed proteins. The absorbance profiles at 280 nm (upper set of peaks) and fluorescence profiles at 475 nm (lower set of peaks) for each elution are shown in Panels A & C. Peak 1 is the nonadsorbed material, and Peak 2 is the eluted target protein. The SDS-PAGE analysis results (using a 4-20% gradient gel for each purification) are shown in Panels B & D. Lane 1: Starting BL21 *E. coli* extract, Lane 2: Nonadsorbed material, Lane 3: Eluted 6xHN-tagged AcGFP1, Lane 4: Starting BL21 *E. coli* extract, Lane M: Molecular weight marker.


그림 4. Use His60 Ni Superflow for the highest binding capacity. The CEP His60 Ni Superflow Starter Kit was used to generate and purify 6xHN-AcGFP1. Lane M: Molecular weight marker, Lane OS: Clarified lysate with expressed 6xHN-AcGFP1, Lane FT: Nonadsorbed material, Lanes W1 and W2: Wash, Lane E1-E4: Eluted 6xHN-AcGFP1 fractions.

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Bacterial Expression and Purification Starter Kits

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
(Superior yields and better purity were obtained when using His60 Ni Superflow for purification as compared to Competitor Q's Ni resin.)
Clontech's pEcoli Linear Expression System was used to express 6xHN-AcGFP1 in *E. coli*. Equivalent amounts of the same sample were used for comparing His60 Ni Superflow Resin purification (Panel A) with Competitor Q Ni Resin purification (Panel B). Lane M: Molecular weight marker, Lane OS: Original sample, Lane FT: Flowthrough, Lanes W1 and W2: Wash, Lane E1-E4: Eluted 6xHN-AcGFP1 fractions.

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