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표 1. Cells transfected with Xfect have High Viability

Cell Line	Xfect (%)	제품 L (%)
HeLa	79.4 ± 17.9	53.2 ± 26.7
HEK-293	63.5 ± 12.5	52.3 ± 12.4
CHO-K1	86 ± 0.8	90.9 ± 10.7
MCF7	51.4 ± 10.6	29.7 ± 21.3


그림 2. Obtain high transfection efficiencies in many cell types with Xfect, including hard-to-transfect cells. Xfect and Competitor Product L (Panel A) or Xfect and Competitor Product LX (Panel B) were used according to their respective protocols to transfect the cell lines listed above with plasmid DNA encoding AcGFP1, in a 6-well format, 48 hr posttransfection, AcGFP1 expression was assessed by flow cytometry in order to determine transfection efficiency.

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
그림 3. Xfect Stem provides superior transfection efficiency in mES cells. Xfect Stem and three competitor reagents were used according to their respective protocols to transfect the cell lines listed above using plasmid DNA encoding AcGFP1, in a 6-well format, 48 hr posttransfection, AcGFP1 expression was assessed by flow cytometry in order to determine transfection efficiency.

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
그림 4. Xfect Stem exhibits far higher transfection efficiency than leading competitor reagents. ES-E14TG2a mES cells were transfected in a 6-well plate format with a plasmid expressing AcGFP1 using Xfect Stem, Product L, Product LX, or Product F according to each manufacturer's recommended protocol. 48 hr posttransfection, the cells were imaged using white light (lefthand images) and by fluorescence microscopy using a Zeiss Axioskop microscope equipped with a GFP filter (righthand images),ansfection efficiency.

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제품명	용량	TaKaRa Code
Xfect	100 rxns	631317
	300 rxns	631318
Xfect Stem	100 rxns	631320
	300 rxns	631321