혁신적인 directional cloning 기술

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
그림 1. The In-Fusion Advantage Cloning Protocol. A thermostable, high-fidelity DNA polymerase and gene-specific primers are used to generate a PCR product containing a gene of interest with 15 bp extensions homologous to the ends of vector. The PCR product is either treated with Cloning Enhancer or spin column-purified, and combined with the vector in a 30 min In-Fusion Advantage cloning reaction. The cloning reaction is then used to transform *E. coli*.

혁신적인 directional cloning 기술

In-Fusion™ Advantage PCR Cloning Kit

Continued...

표 1. In-Fusion PCR Cloning vs. General Cloning Systems

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표 2. Transformation 효율 (Number of Colonies)¹⁾

Insert Size (kb)	In-Fusion ²⁾	A사 S Kit ²⁾	B사 T Kit
None (control)	5	2	9
0.5	~950	19	236
1.0	~1,600	11	164
2.0	~1,570	20	~540
4.0	~1,050	18	298
6.0	~1,480	57	101
8.0	80	0	114
12.0	360	0	84

¹⁾ Cloning was performed according to the instructions for each kit, with transformation into the required E. coli cell lines. The numbers of colonies generated for the indicated inserts are shown. Prior to setting up the cloning reactions, PCR fragments between 0.5 kb and 6 kb in size were spin column-purified, and the 8 kb and 12 kb fragments were gel-purified. Because the Cloning Enhancer was used to treat the 0.5 kb 6 kb PCR inserts for use with the In-Fusion PCR Cloning Kit, there was no need to purify these inserts. All of the transformants were plated on LB/Amp plates 9, 10).

²⁾ For the In-Fusion and S kit products, only 1/10 of the total transformation mixes were plated for the 0.5 kb, 1 kb, 2, kb, 4 kb, and 6 kb insert sizes. For all other fragment sizes and the Directional T kit, the entire transformation was plated.

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혁신적인 directional cloning 기술 In-Fusion™ Advantage PCR Cloning Kit

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표 3. Primers for PCR Amplification of the PD-L2 Extracellular Domain, IgG2a Fc, and IgA Tailpiece

PD-L2 Extracellular Domain with 15 bp Overlaps to NcoI-Digested Vector and IgG2a Fc, 693 bp PCR Product	
Sense (NcoI site is underlined)	5'-TTCAAATCCACCATGGAGACAGACACTCTGCTAT-3'
Antisense	5'-CTTGATTGTGGGCCCTCTGGGGACTTTGGTTCCATCCGA-3'
IgG2a Fc, 678 bp PCR Product	
Sense	5'-GGGCCCAACAATCAAGCCCTGTCTC-3'
Antisense	5'-CCGGGAGAAGCTCTTAGTC-3'
IgA Tailpiece with 15 bp Overlaps to IgG2a Fc and Sall-Digested Vector, 100 bp PCR Product	
Sense	5'-AAGAGCTTCTCCCGGCTGTCTGGGTAACCCACCAATGTCAG-3'
Antisense (Sall site is underlined)	5'-AGTAACGTTAGTCGACTCAGTAGCAGATGCCATCTC-3'


그림 2. Seamless construction of a three-domain immunoglobulin fusion protein with a four-piece In-Fusion reaction. Segments were generated by PCR using primers that contained a 15 bp overlap with the adjacent segment, and 20-30 bp of segment-specific sequence. Colored regions indicate overlap regions with 15 bp of identity. Arrows indicate PCR primers. The segments were joined to a Nco I -Sal I -digested expression vector in a ligase-free In-Fusion reaction.

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혁신적인 directional cloning 기술

In-Fusion™ Advantage PCR Cloning Kit

Continued...

표 4. Transformation Efficiency of a Four-Piece In-Fusion Reaction

PD-L2 Extracellular Domain	IgG2a Fc	IgA Tailpiece	Vector	Total Transformants	No. Error-Free/No. Sequenced
693 bp 6.1 ng	678 bp, 6.0 ng	100 bp 0.9 ng	11,365 bp 50 ng	120	2/2


그림 3. Recombinant Antibody Protein Production.
One μ g of purified protein was analyzed using SDS-PAGE under reducing conditions and stained with Coomassie blue.
Lane 1: Protein molecular weight standards with size in kDa indicated.
Lane 2: Murine PD-L2-IgG2a Fc-IgA tailpiece fusion protein.

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제품명	형태	용량	Cloning Enhancer Included	Spin Columns Included	TaKaRa Code
In-Fusion Advantage PCR Cloning Kit w/Stellar Competent Cells	Liquid	10 rxns			639629
In-Fusion Advantage PCR Cloning Kit w/Cloning Enhancer & Stellar Competent Cells	Liquid	10 rxns	Yes		639630
In-Fusion Advantage PCR Cloning Kit w/NucleoSpin & Stellar Competent Cells	Liquid	10 rxns		Yes	639631
In-Fusion Advantage PCR Cloning Kit	Liquid	10 rxns			639619
	Liquid	50 rxns			639620
	Liquid	100 rxns			639621
In-Fusion Advantage PCR Cloning Kit w/Cloning Enhancer	Liquid	10 rxns	Yes		639616
	Liquid	50 rxns	Yes		639617
In-Fusion Advantage PCR Cloning Kit w/NucleoSpin	Liquid	100 rxns	Yes		639618
	Liquid	10 rxns		Yes	639622
In-Fusion Advantage PCR Cloning Kit w/NucleoSpin	Liquid	50 rxns		Yes	639623
	Liquid	100 rxns		Yes	639624

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- Evans, D.H, *et al.* (2003) In patent application 20030162265, USA, pp. 1- 52.
- Hamilton, M.D, *et al.* (2007) *Nucleic Acids Res*, **35**(1):143 - 151.
- Marsischky, G, & LaBaer, J, (2004) *Genome Res*, **14**(10B):2020-2028.
- In-Fusion™ PCR Universal Cloning Kits, (October 2005) *Clontechniques XX* (2):16- 17.
- In-Fusion™ 2.0 CF Liquid PCR Cloning Kits, (April 2008) *Clontechniques XXIII*(2):20-21.
- Fluorescent Protein Vectors, (April 2008) *Clontechniques XXIII* (2):11-14.
- In-Fusion™ Ready BacPAK Vector Set, (July 2006) *Clontechniques XXI* (2):16-17.
- Benoit, R, M, *et al.* (July 2006) *Clontechniques XXI* (2):21- 22.
- Superior One-Step Cloning of PCR Fragments into Any Vector with the In-Fusion 2.0 Cloning System, (April 2007) *Clontechniques XXII* (2):13-15.
- Efficient Cloning of Long PCR Inserts with the In-Fusion PCR Cloning System, (April 2007) *Clontechniques XXII* (2):16-17.
- Marsischky, G, & LaBaer, J, (2004) *Genome Res*, **14**(10B):2020-2028.
- Hartman, S, *et al.* (January 2005) *Clontechniques XX*(1):26-27.
- Berrow, N,S, *et al.* (2007) *Nucleic Acids Res*, **35** (6):e45.
- Zhu, B, *et al.* (2007) *BioTechniques* **43**(3):354-359.
- In-Fusion™ PCR Cloning Kits-FAQs (July 2007) *Clontechniques XXII* (3):22- 24.
- Design Genes with Ease Using In-Fusion™ Cloning, (April 2008) *Clontechniques XXIII* (2):22-24.
- Latchman, Y, *et al.* (2001) *Nature Immunol*, **2**(3):261-268.
- Hirano, N, *et al.* (2006) *Blood*, **107**(4):1528 - 1536.
- Arthos, J, *et al.* (2002) *J. Biol. Chem*, **277**(13):11456-11464.

* License Notice : [1]