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Molecular Biology Group
Clontech Laboratories, Inc.

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
그림 1. Organization of the anti-lamin A/C and anti-Luc (L1) miRNA expression cassettes, Anti-luciferase and anti-lamin A/C miRNA-like sequences, based on the miR-30 miRNA, were joined to the coding regions of AcGFP1 or DsRed2 cDNAs in the pTRE-Tight expression plasmid. Sense (S) and anti-sense (AS) target sequences are connected by a 19 nt loop sequence. Transcripts are processed to yield either an active miRNA or a mature mRNA for the fluorescent protein.

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
그림 2. Doxycycline-Induced suppression of lamin A/C in HeLa cell lines. Stable, inducible HeLa cell lines were produced that expressed the Tet-On Advanced transactivator and either the pTRE-Tight-lamin A/C-AcGFP1 or DsRed2 construct. Cells were grown in the presence or absence of Dox and assayed after 72 hr by Western Blot analysis. Lane 1: Parental cells. Lane 2: Anti-lamin A/C-AcGFP1 cells without Dox. Lane 3: Anti-lamin A/C-AcGFP1 with Dox (1 µg/ml). Lane 4: Anti-lamin A/C-DsRed2 without Dox. Lane 5: Anti-lamin A/C-DsRed2 with Dox (1 µg/ml).

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
그림 3. Doxycycline-induced fluorescence in anti-lamin A/C HeLa cell lines. The HeLa cell lines described and treated as in Figure 2 were visualized by fluorescence microscopy. Panel A, Anti-lamin A/C-AcGFP1 cells without Dox, Panel B, Anti-lamin A/C-AcGFP1 cells with Dox (1 ug/ml), Panel C, Anti-lamin A/C-DsRed2 cells without Dox, Panel D, Anti-lamin A/C-DsRed2 cells with Dox (1 ug/ml). Dox treatment-induced strong fluorescence in cells that also exhibited lamin A/C knockdown.


그림 4. Doxycycline treatment reduces luciferase activity in HEK-293 cells transfected with the Anti-Luc-DsRed2 reporter. Hek-293 cell lines stable transformed with Tet-On Advanced were transiently cotransfected with pCMV-Luc and pTRE-Tight-L1-DsRed2 in triplicate at a ratio of 1:2, and then grown in the presence of Dox (1 ug/ml) for 72 hr. The L1 miRNA induced by Dox treatment caused a 85% reduction in luciferase activity, S.D. = standard deviation.

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제 품 명	용 량	TaKaRa Code
Tet-On Advanced Inducible Gene Expression System	each	630930
HEK 293 Tet-On Advanced Cell Line	each	630931
HepG2 Tet-On Advanced Cell Line	each	630932
pTRE-Cycle1	20 µg	631061
pTRE-Cycle2	20 µg	631066
pTRE-Cycle3	20 µg	631064

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

Living Colors[®] Antibody

- Western Blot
- Immuno precipitation
- Immuno cytochemistry

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Western blot detection of DsRed and its variants.

A HeLa cells were transiently transfected with the following vectors, and then analyzed by Western blot with the DsRed Monoclonal Antibody. Lane 1: pDsRed2-N1, Lane 2: pDsRed2-C1, Lane 3: pDsRed1-C1, Lane 4: pDsRed2-Mito (N-terminal fusion), Lane 5: pDsRed2-Nuc (C-terminal fusion), Lane 6: control (untransfected cells). A 20 kDa, non-specific band appears with extended exposure.

B Western blot with the DsRed Polyclonal Antibody. Lane 1: 25 ng of recombinant AcGFP1 protein. Lanes 2-5: HEK-293 cells were stably transfected with the following vectors, Lane 2: pAcGFP1, Lane 3: Control (untransfected cells), Lane 4: pDsRed-Monomer-N1, Lane 5: pDsRed-Express-N1.

형광단백질	추천하는 antibody	TaKaRa Code	적용가능한 실험
AcGFP1	GFP Monoclonal Antibody	632375	WB, IP, IC
	A.v. Monoclonal Antibody (JL-8)	632380, 632381	WB, IP, IC
	Full-Length A.v. Polyclonal Antibody	632459, 632460	WB, IP
	A.v. Peptide Antibody	632376, 632377	WB, IP, IC
AmCyan1	Anti-RCFP Polyclonal Pan Antibody	632475	WB
AsRed2	Anti-RCFP Polyclonal Pan Antibody	632475	WB
mBanana	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
mCherry	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
DsRed2	Anti-RCFP Polyclonal Pan Antibody	632475	WB
	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
DsRed-Express	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
DsRed-Express2	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
DsRed-Monomer	DsRed Polyclonal Antibody	632496	WB, IP, IC
HcRed1	Anti-RCFP Polyclonal Pan Antibody	632475	WB
	HcRed Polyclonal Antibody	632452	WB
mOrange	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
mPlum	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
mRaspberry	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
mStrawberry	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
tdTomato	DsRed Monoclonal Antibody	632392, 632393	WB, IP
	DsRed Polyclonal Antibody	632496	WB, IP, IC
ZsGreen1	Full-Length ZsGreen Polyclonal Antibody	632474	WB
	Anti-RCFP Polyclonal Pan Antibody	632475	WB
ZsYellow1	Anti-RCFP Polyclonal Pan Antibody	632475	WB
EGFP*	GFP Monoclonal Antibody	632375	WB, IP, IC
EYFP*	A.v. Monoclonal Antibody (JL-8)	632380, 632381	WB, IP, IC
EYFP*	Full-Length A.v. Polyclonal Antibody	632459, 632460	WB, IP

WB= Western Blot. IP= Immunoprecipitation IC= Immunocytochemistry * 종매품 : vector만 종매