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- Tight Transcriptional Regulation - Tetracycline (Tet)-inducible Expression
- Rapid, reversible Post-translational Regulation-ProteoTuner-controlled Protein degradation
- Control the Level of TWO protein - Independently

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
그림 3. Regulating protein levels is easy with pTRE-Cycle Vectors. HEK 293 cells were cotransfected with pTet-On Advanced and pTRE-Cycle2-ZsGreen1 (encoding mCherry and DD-tagged ZsGreen1). Transfected cells were subsequently grown in medium: -Dox/-Shield1 (Panel A); +Dox/-Shield1 (Panel B); and +Dox/+Shield (Panel C). Cells incubated in medium containing only Dox accumulated only mCherry, whereas cells incubated in medium containing both Dox and Shield1 accumulated both mCherry and DD-tagged ZsGreen1.

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
그림 4. Get tight, independent control over levels of coexpressed proteins. HEK 293 cells were cotransfected with pTet-On Advanced and pTRE-Cycle1-Sox2-Oct4 (encoding both Sox2 and DD-tagged Oct4), and grown on medium containing various concentrations of Dox and Shield1. 48 hrs post-transfection, cell lysates were subjected to western analysis with an anti-Sox2 monoclonal antibody (mAb) and an anti-DD mAb. Cells incubated in medium containing only Dox, accumulated only Sox2 in a concentration-dependent manner (Lanes 2&3); cells incubated in medium containing both Dox and Shield1 accumulated both Sox2 and DD-tagged Oct4, also in a concentration-dependent manner (Lanes 4&5).


그림 5. Protein levels can be adjusted in less than 4hrs. HEK 293 cells were cotransfected with pTet-On Advanced and pTRE-Cycle2-ZsGreen1 (encoding both mCherry and DD-tagged ZsGreen1), and incubated overnight in medium containing Dox and Shield1. Transfectants were subsequently incubated in medium containing only Dox (no Shield1). After 4hrs, mCherry levels were unaffected (Panel A); however, DD-tagged ZsGreen1 was largely degraded (Panel B).

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