Low background, High Sensitivity의 새로운 Inducible System 출시!

# Tet-On<sup>®</sup> 3G Inducible Expression System

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
그림 1. The Tet-On 3G systems allow inducible gene expression only in the presence of doxycycline. When Dox binds, the transactivator undergoes a conformational change allowing it to bind tet operator (*tetO*) repeats within the TREG Promoter ( $P_{TRE3G}$ ). The transactivator activates expression through transcription activation domain repeats.

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
그림 2. The  $P_{TRE3G}$  promoter results in significantly reduced basal expression, HEK 293 cells were transiently cotransfected with both the response vectors (containing luciferase) and regulator vectors from each of the Tet-On 3G and Tet-On Advanced Inducible Expression Systems. The cells were cultured in the presence and absence of Dox, and after 24 hr, luciferase expression was measured. Although both systems provided strong expression in the presence of Dox, the Tet-On 3G System produces far lower background expression in the absence of Dox (inset).


그림 3. Tet-On 3G demonstrates higher sensitivity to doxycycline than Tet-On Advanced. Tet-On 3G and Tet-On Advanced genes were integrated at the same locus in a stable HLF33 cell line expressing luciferase from a TRE promoter. For each of these two double-stable cell lines, induced luciferase expression was measured in response to a range of doxycycline (Dox) concentrations. At 5–10 ng/ml Dox, induced expression was 100–150-fold higher for the Tet-On 3G cell line, and at 50 ng/ml, expression was 4.6 fold higher (data kindly provided by Professor W. Hillen and Dr. C. Berens, University of Erlangen).


그림 4. Tet-On 3G is highly sensitive to as little as 10 ng/ml of doxycycline. Following cotransient transfection of pCMV-Tet3G and pTRE3G-Luc in HeLa cells, increasing levels of Dox were added and expression of luciferase was measured using an anti-luciferase antibody (Panel A) and a luciferase assay (Panel B). Induced expression was very high even with Dox concentrations as low as 10 ng/ml, and was detectable even at 0.1 ng/ml Dox.

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표 1. Three Generations of Tet-On Inducible Expression Systems

Name	Generation	Transactivator Protein	Inducible Promoter
Tet-On System	1st	Tet-On	P <sub>TRE2</sub>
Tet-On Advanced System	2nd	Tet-On Advanced	P <sub>TIGHT</sub>
Tet-On 3G System	3rd	Tet-On 3G	P <sub>TRE3G</sub>

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
그림 5. Sequence comparison of the three generations of TRE promoters sold by Clontech. Each promoter consists of 7 identical repeats of the tet operator sequence (green), and contains a minimal CMV promoter-although this sequence has been mutated to consensus in P<sub>TRE3G</sub><sup>3)</sup>. Moreover, in P<sub>TRE3G</sub> the spacer sequences between the tetO repeats are evenly spaced and contain randomized central sequences.


그림 6. Comparison of the three generations of Tet-On. Panel A shows the transactivator and inducible promoter combinations for each generation of the Tet-On system. Tet-On was launched by Clontech in 1996 at the time this was the premier inducible expression system, and its performance has only been surpassed by subsequent generations of the Tet-On system. Panel B shows co-transient transfection experiments with both vectors in HEK 293 and HeLa cells. Although Tet-On Advanced shows great improvement over the original Tet-On System when comparing the difference between the induced and the uninduced states (fold induction), Tet-On 3G shows even higher fold induction due to the significantly reduced basal expression provided by the P<sub>TRE3G</sub> promoter.

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
그림 10. Tet-On 3G (EF1 $\alpha$  Version) provides inducible expression in hematopoietic cells. Jurkat cells were transfected with pEF1 $\alpha$ -Tet3G using Xfect transfection reagent, and stable clones were selected by limiting dilution. 18 stable clones were then tested for inducibility by transient transfection using pTRE3G-expressing luciferase. Six out of eighteen clones showed more than 2,000 fold induction via transient transfection.


그림 11. Tet-On 3G (EF1 $\alpha$  Version) provides inducible expression in stem cells. Mouse embryonic stem cells (ES-E14TG2a mES cells) were cotransfected with pTRE3G-ZsGreen1 and pEF1 $\alpha$ -Tet3G using Xfect Stem transfection reagent. The stem cells show ZsGreen1 expression only in the presence of Dox (Panel C).

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
그림 12. Co-inducible expression of two genes from a single bicistronic transcript. pTRE3G-IRES expressing both mCherry, from MCS I, and luciferase, from MCS II, was cotransfected with pCMV-Tet3G into HeLa cells with increasing levels of Dox. Induced expression of mCherry was monitored by fluorescent microscopy (Panel A) and luciferase expression was measured using a luciferase assay (Panel B).


그림 13. Monitor inducible expression and easily screen for inducible clones using bright fluorescent proteins. Clones created using the Tet-On 3G Inducible Expression System (with mCherry) or Tet-On 3G Inducible Expression System (with ZsGreen1) will respectively fluoresce red and green, but only in the presence of Dox. The size (bp) of the fluorescent protein gene is optimal for maximum expression of your GOI from the internal ribosome entry site (IRES).

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제품명	용량	TaKaRa Code
Tet-On 3G Inducible Expression System	each	631168
Tet-On 3G Inducible Expression System (EF1 $\alpha$ Version)	each	631167
Tet-On 3G Inducible Expression System (Bicistronic Version)	each	631166
Tet-On 3G Inducible Expression System (with mCherry)	each	631165
Tet-On 3G Inducible Expression System (with ZsGreen1)	each	631164
Tet System Approved FBS	50 ml	631107
	500 ml	631106
Doxycycline	5 g	631311
TetR Monoclonal Antibody (Clone 9G9)	40 $\mu$ g	631131
	200 $\mu$ g	631132

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그림 14. To ensure low background with Tet-On 3G, it is essential to use fetal bovine serum that is functionally tested. Clontech offers four such serum options, including serum that is sourced from the US, Mexico, and Australia, as well as US-sourced serum that is additionally tested for use in mouse embryonic stem cells.

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- Low, R., Heinz, N., Hampf, M., Bujard, H. & Gossen, M. Ameliorating the dynamic properties of the Tet system by altered minimal promoter design. *Manuscript submitted for publication.*
- Zhou, X., Vink, M., Klaver, B., Berkhout, B. & Das, A. T. Optimization of the Tet-On system for regulated gene expression through viral evolution. (2006) *Gene Ther.* **13**(19):1382 – 1390.
- Gossen, M. & Bujard, H. (1992) *Proc. Natl. Acad. Sci. USA* **89**(12):5547 – 5551.
- Urlinger, S. *et al.*, (2000) *Proc. Natl. Acad. Sci. USA* **97**(14):7963 – 7968.

\*License Notice : [3]