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
그림 1. Typical gel profile of ds cDNA synthesized with the SMARTer™ PCR cDNA Synthesis Kit (Panel A) and the SMARTer™ Pico PCR cDNA Synthesis Kit (Panel B), using the Control Human Placental Total RNA as a template. 2 ng (Panel A) or 1 ng (Panel B) of the Control Human Placental Total RNA was subjected to first-strand cDNA synthesis and purification as described in the user manual for each kit. 10 μl was used for SMARTer™ PCR amplification (Panel A) and 80 μl was used for SMARTer™ Pico PCR amplification (Panel B). A range of PCR cycles was performed (18, 21, 24 and 27). 5 μl of each PCR product was electrophoresed on a 1.2% agarose/EtBr gel in 1X TAE buffer following the indicated number of PCR cycles. Lanes M : 1 kb DNA ladder size markers. The arrows indicate the strong band at 900 bp typically seen for human placental total RNA.

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SMARTer	SMARTer Pico
<ul style="list-style-type: none"> <li>• 2~1000ng total RNA</li> <li>• Template reaction volume up to 3.5 μl</li> <li>• Total RNA template concentration as low as 0.6 ng/μl</li> </ul>	<ul style="list-style-type: none"> <li>• 1~1000ng total RNA</li> <li>• Template reaction volume up to 50 μl</li> <li>• Total RNA template concentration as low as 20 pg/μl</li> </ul>
<ul style="list-style-type: none"> <li>• SMARTer first-strand cDNA synthesis</li> <li>• Volume=10 μl</li> <li>• Dilute 1:5 with TE Buffer</li> <li>• Volume=50 μl</li> </ul>	<ul style="list-style-type: none"> <li>• SMARTer Pico first-strand cDNA synthesis</li> <li>• Volume=106 μl</li> <li>• Purify with NucleoSpin® Column</li> <li>• Column Elution Volume = 80 μl</li> </ul>
<ul style="list-style-type: none"> <li>• Use 10 μl cDNA for SMARTer PCR amplification</li> <li>• 100 μl reaction</li> <li>• Cycle optimization and scale-up</li> </ul>	<ul style="list-style-type: none"> <li>• Use 80 μl cDNA for SMARTer Pico PCR amplification</li> <li>• 100 μl reaction</li> <li>• Cycle optimization and scale-up</li> </ul>
<ul style="list-style-type: none"> <li>• Purify PCR products with NucleoSpin®</li> <li>• Yields 1-2 μg ds cDNA</li> </ul>	<ul style="list-style-type: none"> <li>• Purify PCR products with NucleoSpin®</li> <li>• Yields 1-2 μg ds cDNA</li> </ul>

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
그림 2. The SMARTer II A Oligo increases the overall yield of the 5' -RACE PCR products. Note the different intensities of the 2.6 kb 5' -RACE product in Lane 2 (produced by the SMART II A oligo) versus that in Lane 6 (produced by the SMARTer II A oligo); both bands resulted from the amplification of only 2 ng of original RNA template. In the assay shown, 2 ng and 50 ng of human placental total RNA, 12 pmol SMART II A oligo or SMARTer II A oligo, and 100 units of SMARTScribe Reverse Transcriptase were combined in a first-strand cDNA synthesis reaction. The resulting single-stranded cDNA was then used as a template for 5' -RACE PCR amplification according to the protocol in the user manual (PT4096-1). The 5' -RACE PCR reaction was performed using the Universal Primer Mix and a gene-specific primer designed for transferrin receptor (TFR). Lanes 1 & 5: No RNA 5' -RACE control. Lanes 2 & 6: TFR 5' -RACE from cDNA synthesized using 2 ng of total RNA as template. Lanes 3 & 7: TFR 5' -RACE from cDNA synthesized using 50 ng of total RNA as template. Lanes 4 & 8: 5' -cDNA internal control using 2 ng total RNA as template. Lane M: 1 kb DNA Ladder.

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1. Zhu, et al. (2001) *BioTechniques* 30:892 - 897.
2. Chenchik, A. et al. (1998) In *Gene Cloning and Analysis by RT-PCR*. Eds Siebert, P. & Larrick, J. (BioTechniques Books, MA) Ch 22.
3. Be SMART<sup>™</sup> About First-Strand cDNA Synthesis (January 2009) *Clontechiques* XXIV(1):15-17.

\*License Notice : [6]


그림 3. SMARTer 5' - and 3' -RACE results for different genes. In this assay, we used the SMARTer RACE Kit to amplify 5' - and 3' -RACE products from many different genes, starting with human placental total RNA. The gel shows several representative 5' - and 3' -RACE products amplified from cDNA encoding the following: Lane 1: 1 kb DNA marker, Lane 2: 100 bp DNA marker, Lanes 3 & 4: tyrosyl-tRNA synthetase (YARS), Lanes 5 & 6: interferon gamma receptor 1 (IFNGR1), Lanes 7 & 8: adrenergic, beta, receptor kinase 1 (ADRBK1), Lanes 9 & 10: actin, beta (ACTB), Lanes 11 & 12: transferrin receptor (TFR), Lane 13: low density lipoprotein-related protein 2 (LRP2). The 1,3 kb 5' -RACE PCR product shown was amplified from cDNA synthesized from the gene's 15.7 kb mRNA. Lane 14: low density lipoprotein-related protein 1 (LRP1). The 2.1 kb 5' -RACE PCR product shown was amplified from cDNA synthesized from the gene's 14.9 kb mRNA.

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