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
그림 2. High Performance Plasmid Transfection. Primary Human Small Epithelial Cells (HSAEpic) were transfected using TransIT-2020 and an EGFP expression plasmid (4:1 reagent-to-DNA ratio). Images were taken 24 hours post-transfection using a Zeiss axiovert inverted fluorescence microscope.

■ TransIT[®]-2020 시리즈

제품명	용량	TaKaRa Code
TransIT [®] -2020 Transfection Reagent	0.4 ml	V5404
TransIT [®] -2020 Transfection Reagent	1 ml	V5400
TransIT [®] -2020 Transfection Reagent	5 x 1 ml	V5405
TransIT [®] -2020 Transfection Reagent	10 x 1 ml	V5406

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Experimental Details: Luciferase gene expression was compared in HEK 293, K562, CHO-K1, and primary MEF cells transfected with a luciferase expression plasmid using TransIT-2020 (Mirus Bio, 3:1 reagent to DNA ratio), FuGENE HD (Roche, 7:2 reagent to DNA ratio), Lipofectamine[™] 2000 (Invitrogen, 5:2 reagent to DNA ratio), and Lipofectamine[™] 2000 CD (Invitrogen, 5:2 reagent to DNA ratio). Transfections were performed in 24-well plates using 0.5 μ g of plasmid DNA per well. The optimal level of each transfection reagent was determined empirically and all reagents were used according to manufacturer's protocol. Cells were harvested at 24 hours post-transfection and assayed for luciferase activity, Error bars represent the standard deviation of triplicate wells.

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Cell Line	Description	사용 연구자 검증	Mirus 검증
2B4.11	Mouse T-lymphocyte hybridoma cell line	○	
A549	Human lung carcinoma cell line	○	
C3H/10T1/2	Clonal mouse embryo cell line	○	
C6	Rat glioma cell line		○
CCD-18Co	Human colonic myofibroblast cell line	○	
CHO-K1	Chinese hamster ovary cell line	○	○
Colon 38	Mouse colon adenocarcinoma cell line	○	
COS-7	Monkey kidney fibroblast-like cell line	○	○
Daoy	Human medulloblastoma cell line		○
DB-TRG-05MG	Human glioma cell line		○
DI-TNC1	Rat astrocytoma cell line		○
HCT 116	Human colorectal carcinoma cell line	○	
HEK-293	Human embryonic kidney cell line	○	○
HEK-293F	Human embryonic kidney cell line	○	
HEK-293T	Human embryonic kidney cell line	○	
HeLa	Human cervix adenocarcinoma		○
HepG2	Human hepatocellular carcinoma cell line		○
HL-1	Mouse atrial cardiomyocyte tumor cell line	○	
HPTC	Human thyroid carcinoma cell line	○	
INS-1	Rat insuloma cell line	○	
J774A.1	Mouse monocyte/macrophage sarcoma cell line	○	
Jurkat	Human T-lymphocyte leukemia cell line		○
K-562	Human myelogenous leukemia cell line		○
KYSE 410	Human oesophageal cancer cell line	○	
LS180	Human colon adenocarcinoma cell line	○	
MCF-7	Human mammary gland adenocarcinoma cell line	○	
MDA-MB-231	Human mammary gland adenocarcinoma cell line	○	
Monc-1	Mouse neural crest stem cell line	○	
NCCIT	Human teratocarcinoma cell line	○	
NCI-H522	Human lung adenocarcinoma cell line	○	
NE-1	Human oesophageal cancer cell line	○	
Neuro-2a	Mouse neuroblastoma cell line		○
NIH-3T3	Mouse embryonic fibroblast cell line	○	○
NRK-52E	Rat kidney epithelial cell line	○	
PC-12	Rat adrenal gland pheochromocytoma cell line	○	
Phoenix [™] Eco	Human embryonic kidney cell line	○	
RAW 264.7	Mouse monocyte/macrophage tumor cell line	○	
RGC-5	Rat retinal ganglion cell line	○	
Saos-2	Human osteosarcoma cell line	○	
SH-EP1	Human neuroblastoma cell line	○	
SH-SY5Y	Human neuroblastoma cell line	○	
SK-N-BE(2)	Human neuroblastoma cell line	○	
T-98G	Human glioblastoma cell line	○	
TE-8	Human oesophageal cancer cell line	○	
U2OS	Human osteosarcoma cell line	○	
Vero	Monkey kidney epithelial cell line		○

■ **TransIT[®]-2020 transfection 시약이 적용된 primary와 stem cell line**

Primary or Stem Cells	사용 연구자 검증	Mirus 검증
Human airway epithelial cells	○	○
Human astrocytes	○	
Human ES cell derived neural progenitors	○	
Human mesenchymal stem cells	○	
Human osteoblasts	○	
Human skin fibroblasts	○	
Mouse cortical neurons (E18)	○	
Mouse embryo fibroblasts	○	○
Mouse ES cell derived cardiomyocytes	○	
Mouse osteoblasts	○	
Rat aortic smooth muscle cells (SV40 Large T-antigen transformed)	○	