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# SMARTer<sup>®</sup> Universal Low Input RNA Kit for Sequencing

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

그림 2. Synthesis scheme for SMARTer Universal Low Input RNA Kit for Sequencing. First-strand cDNA is synthesized from sheared rRNA depleted RNA, using a modified N6 primer called SMART N6 CDS (N = A, G, T, or C). When SMARTScribe Reverse Transcriptase reaches the 5' end of the RNA, its terminal transferase activity adds a few additional nucleotides to the 3' end of the cDNA. This non-template nucleotide stretch anneals with the SMARTer Oligonucleotide, creating an extended template where the RT continues through the II A sequence. The SMARTer anchor sequence and the modified N6 sequence serve as universal priming sites for cDNA amplification. Finally, amplified cDNA is digested with RsaI to remove the adapter prior to sequencing.

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**그림 3. Yield and purity of input RNA and resulting cDNA libraries.** Panel A. Human Brain Total RNA (HBTR) was chemically sheared, spiked with ERCC control RNA (4 µl of a 1:1,000 dilution per 100 ng), and rRNA-depleted using a modified Ribo-Zero protocol for lowinput samples. One microliter was analyzed using the Agilent 2100 Bioanalyzer (RNA 6000 Pico chip). Panel B. Three amounts of rRNAdepleted RNA (156 pg, 780 pg, and 7.8 ng) were used as input for the SMARTer Universal Low Input RNA Kit for Sequencing, and 1 µl of the resulting cDNA library was analyzed on the Agilent 2100 Bioanalyzer (High Sensitivity DNA chip).

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| Table 1: Sequence Alignment Metrics |                  |                  |                  |
|-------------------------------------|------------------|------------------|------------------|
| Input RNA                           | 156 pg           | 780 pg           | 7.8 ng           |
| No. of Reads                        | 27,043,029       | 25,247,363       | 16,991,089       |
| Mapped to rRNA                      | 316,377 (1.2%)   | 360,696 (1.4%)   | 324,628 (1.9%)   |
| Mapped to Mitochondrial RNA         | 6,408,693 (24%)  | 6,007,327 (24%)  | 3,858,048 (23%)  |
| Mapped to RefSeq                    | 16,863,730 (84%) | 16,158,298 (86%) | 10,977,858 (87%) |
| Uniquely Mapped to RefSeq           | 15,480,061 (77%) | 14,823,088 (79%) | 10,034,740 (79%) |
| Exons                               | 6,464,355 (38%)  | 6,185,920 (38%)  | 4,256,916 (39%)  |
| Introns                             | 10,399,375 (62%) | 9,972,378 (62%)  | 6,720,942 (61%)  |
| Genes Identified                    | 15,140           | 15,574           | 15,697           |

**표 1. Sequencing alignment metrics for three N6 libraries.** cDNA libraries were created using 156 pg, 780 pg, and 7.8 ng of rRNA-depleted human brain reference RNA and used for sequencing. Illumina adapters and indices were added to 1-2 ng of each cDNA library using the Low Input Library Prep Kit, and the cDNA samples were sequenced on an Illumina MiSeq Platform with 1 x 50 bp reads. Reads were trimmed by CLC Genomics Workbench and mapped to rRNA and the mitochondrial genome with CLC (% reads indicated). The unmapped reads were subsequently mapped with CLC to the human genome with the RefSeq masking, producing mapped reads and uniquely mapped reads. The number of genes identified in each library was determined by the number of genes with an RPKM (read per kilobase of exon per million of reads) of at least 0.1. The number of reads that map to introns or exons is a percentage of the reads successfully mapped to RefSeq. Note: The 156 pg, 780 pg, and 7.8 ng of rRNAdepleted RNA correspond to 2 ng, 10 ng, and 100 ng of total RNA, prior to rRNA depletion.

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
그림 4. High reproducibility and sensitivity across a wide RNA input range. Scatter plots of expression (RPKM, read per kilobase of exon per million of reads) for cDNA libraries prepared from rRNA-depleted human brain RNA. Panels A-C. Comparisons of pairs of cDNA library replicas (Replica #1 and Replica #2) created from 156 pg, 780 pg, and 7.8 ng of input RNA show high reproducibility across a wide range of RNA concentrations. Panels D-F. Comparisons of cDNA libraries generated from pairs of RNA input amounts (780 pg and 156 pg, 7.8 ng and 156 pg, and 7.8 ng and 780 pg) show a high correlation, suggesting consistency across input levels. Axes are plotted on a log<sub>10</sub> scale. Insets indicate the coefficient of correlation by Spearman analysis ( $\rho$ ) and by Pearson correlation ( $R$ ).

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


그림 5. High reproducibility confirmed by ERCC analysis. The ERCC (External RNA Control Consortium) Spike-In Mix (92 transcripts) was chemically sheared and 4 microliters (1:1,000 dilution) was spiked into 100 ng of chemically sheared human brain total RNA prior to rRNA depletion. Two cDNA library replicates were generated with each of 156 pg, 780 pg, and 7.8 ng of rRNA-depleted RNA/ERCC mix using the SMARTer Universal Low Input RNA Kit for Sequencing. RPKM of the ERCC transcripts was plotted against the attomoles per transcript spiked in. The slope, R<sup>2</sup>, and the number of transcripts are indicated per replicate. Pearson correlation ( $R$ ) values per replica set were determined to be 0.998, 0.999, and 0.992 for Panels A-C, respectively. Axes are plotted on a log<sub>2</sub> scale.

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**그림 6. Uniform gene body coverage.** Overlaid data comparing transcript coverage from cDNA libraries generated using 156 pg, 780 pg, and 7.8 ng of rRNA-depleted human brain RNA. The x-axis represents gene length (RefSeq) normalized to 100%, where 0 is the 5'-end of each transcript and 100 is the 3'-end. The y-axis represents the read coverage relative to the highest coverage percentile.

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**그림 7. High correlation with qPCR data via MAQC analysis.** Scatter plots were used to compare differential expression data obtained by sequencing cDNA libraries created with the SMARTer Universal RNA Kit [106 pg HURR/159 pg HBRR (A), 530 pg HURR/790 pg HBRR (B), and 5.3 ng HURR/7.9 ng HBRR (C)] and quantitative PCR (qPCR) data available for HURR and HBRR through the MAQC (MicroArray Quality Control) project (2). The slope and correlation (R) for the comparison line of expression ratio (in RPKM) and qPCR ratio (in Ct) are plotted for HURR and HBRR, in a log<sub>2</sub> scale. The number of genes present in both data sets is indicated. HURR = Human Universal Reference RNA, HBRR= Human Brain Reference RNA.

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