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
그림 1. Entry Points for Transfection. Adult fibroblast cells can be transfected or transduced via several methods (e.g. recombinant virus, plasmid, protein, mRNA, small molecule and miRNA) with a combination of transcription factors including KLF4, SOX2, c-MYC, NANOG, OCT-4 and LIN-28 to reprogram the cells to a pluripotent state. iPS cells can then be differentiated to a myriad of cell types through growth factor addition and/or transfection of selection markers driven by cell type specific promoters. Stem cell derived cell types such as cardiomyocytes, adipocytes, neural cells, pancreatic b-cells, and hematopoietic progenitor cells provide researchers with relevant models for their experiments.

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
그림 2. Fibroblast Transfection with *TransIT*-mRNA. The *TransIT*-mRNA Transfection Kit was used to transfect BJ human neonatal foreskin fibroblasts (A) and MRC-5 human lung fibroblasts (B) with a pseudouridine and 5mC modified based GFP mRNA (Trilink Biotechnologies, Inc.). Transfections were performed in 12-well plates using 1-3 μ l of *TransIT*-mRNA Transfection Reagent and mRNA Boost Reagent to deliver 1 μ g of RNA (1:1:1, 2:2:1 and 3:3:1; reagent: boost: RNA ratio). Cells were assayed 18 hours post-transfection on a Guava HT easyCyte flow cytometer. Toxicity was measured using propidium iodide stain (black line).

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
그림 3. High Efficiency Transfection and Electroporation of Human iPS Cells. The *TransIT*-2020 Transfection Reagent was used to transfect 0.5×10^6 iPS cells with a ZsGreen expression plasmid (Clontech) (A). Transfections were performed in 6-well plates using 7.5 μ l of *TransIT*-2020 Transfection Reagent to deliver 2.5 μ g of DNA (3:1; reagent: DNA). The Ingenio Electroporation Kit was used to transfect 2×10^6 iPS cells on the Amaxa Nucleofector II Device (B). Cells were electroporated with 8 μ g ZsGreen expressing plasmid (Clontech) in 100 μ l and plated in 6-well plates at 0.33×10^6 cells/well. Cells were visualized 24 hours post-transfection and imaged at 4X objective with an Olympus IX71 Inverted Microscope. Images were acquired using phase contrast and green fluorescence. Cells were assayed 24 hours post-transfection on an Accuri Cytometer. The histogram shows untransfected cells (black line) compared to cells transfected with plasmid (green line).

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
그림 4. Plasmid DNA Delivery to iCell Cardiomyocytes using *TransIT-LT1*. Panel A illustrates high efficiency transfection of a GFP encoding plasmid. iCell Cardiomyocytes were plated at 20,000 cells/well in a 96 well tissue culture plate coated with 0.1% gelatin.

After allowing the cells to recover from thaw, cells were transfected with 100 ng/well of pMAXGFP(Amaxa) using *TransIT-LT1* Transfection Reagent with a 2:1 (reagent:DNA) ratio according to the manufacturer's instructions. Fluorescent images were taken 3 days post transfection.

Panel B is a schematic of agonist binding inducing G protein (Gs) mediated activation of adenylyl cyclase which converts ATP to cAMP. The second messenger is able to bind to protein kinase A (PKA) and lead to phosphorylation of the cAMP response element-binding protein (CREB) protein. Upon translocation to the nucleus CREB is able to bind the cAMP response element (CRE) and initiate expression of the luciferase reporter.

Panel C illustrates cAMP induction measured via a luciferase reporter plasmid. iCell Cardiomyocytes were plated for 5 days and subsequently replated using 40,000 or 80,000 cells/well in a 96 well plate pre-coated with gelatin. Three days post-replating cells were transfected using *TransIT-LT1* and the CRE-luciferase reporter plasmid pGL4.29 (Promega). After 18 hours the cAMP pathway was induced using 10 mM isoproterenol for 6 hours. Luciferase activity was measured using the Promega Dual Glo Luciferase Assay. Data is normalized to the control reporter plasmid pGL4.75 (Promega).


그림 5. Efficient siRNA-mediated Gene Silencing by *TransIT-TKO* in iCell Cardiomyocytes. Panels A and B show the effect of GAPDH-targeted siRNA on GAPDH (targeted) and HPRT1 (non-targeted) mRNA expression, respectively. iCell Cardiomyocytes were cultured for 7 days in a 12-well cell culture plate before transfection with either control (scrambled) or GAPDH siRNA (sense: GCUCAUUUCUGGUAUGACUU; antisense: GUCAUACCAGGAAAUGAGCUU) using *TransIT-TKO* (3 – 5 µl/well). 72 hours post-transfection the GAPDH and HPRT1 (non-targeted) mRNA levels were measured relative to 18s rRNA levels and normalized to the mRNA levels obtained following transfection of the control siRNA in each experiment. The bar graphs show the mean with standard error of the mean (SEM) of 3 independent transfection complexes

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
Enables superior gene expression in a variety of cell types. The TransIT-X2™ Dynamic Delivery System and Lipofectamine® 2000 Transfection Reagent were used to transfect plasmid DNA encoding luciferase into 30 different cell types at three reagent-to-DNA ratios. Luciferase expression was compared at 24 hours post-transfection using a standard luciferase assay. Head-to-head comparisons illustrate superior or equal luciferase expression using TransIT-X2 in 26 of 30 cell types; 11 cell types had expression levels 2-fold higher than Lipofectamine 2000 (denoted with †).

† Cell types with >2-fold luciferase expression in head-to-head comparisons.