## Technical Note 1

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


**그림 1.** Size distribution of amplified cDNA across total RNA input. Amplified cDNA using 50, 20, and 5 ng of total HeLa RNA show consistent profile and size distribution across RNA input range upon analysis of Agilent Bioanalyzer traces.

**그림 2.** qPCR and yields reproducibility across total RNA input and transcript abundance. 50, 20, and 5 ng of HeLa total RNA were amplified in triplicate. qPCR results for 7 transcripts show reproducible and consistent results with and without purification of amplified cDNA. Yields, shown on the left axis of the graph, also show highly reproducible degrees of amplification across inputs

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
**그림 3.** qPCR and yields reproducibility across multiple product lots. 5 ng of UHR total RNA was amplified in triplicate for 3 different lots of the Ovation RNA Amplification System V2. qPCR results for 5 genes show reproducible and consistent amplification. Amplified cDNA yields, shown on the left axis of the graph, also show very consistent amplification across multiple lots.

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**그림 4.** qPCR data shows expected and consistent fold change: Using a tissue mixing model, amplification of 5 ng of a 100% Placenta RNA sample vs. a 50% mix of Placental and Spleen RNA and consequent qPCR assays of a panel of 18 Placenta-specific genes showed that the expected 2 fold changes observed and maintained across the different transcripts. The rightmost 4 genes in the graph are housekeeping genes, and the very rightmost gene, RPL35-1 was used as a normalization gene.

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
| Raw Q | Scaling Factor | Back-ground | % Present | (3'/5') GAPDH | 3'/5' Actin |
|-------|----------------|-------------|-----------|---------------|-------------|
| 0.97  | 0.69           | 34.5        | 64.1      | 1.30          | 4.1         |
| 0.98  | 0.59           | 34.4        | 65.6      | 1.21          | 3.8         |

and labeled in duplicate, and hybridized to HG-U133A 2.0 GeneChip arrays, show robust and reproducible array metrics.

**그림 6.** High reproducibility demonstrated by signal correlation and call concordance : 5 ng of HeLa RNA, amplified, fragmented and labeled in duplicate, as previously described show high signal correlations of R<sup>2</sup> of 0.994 and call concordance of 93.2% demonstrating very high level of reproducibility between independently processed samples.

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