Technical Note 2

WT-Ovation™ FFPE System

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표 1. A. Affymetrix HG-U133A 2.0 array performance metrics. Listed are the number of arrays (N), Scaling Factor (SF), Background (Bkgd), and Percent of Present calls (%P) with standard deviation listed in parenthesis (sd) for each value. B shows the array average Pearson correlation matrix, demonstrating the high level of reproducibility between technical replicates.

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
그림 1. Venn diagram of concordance of present calls between FFPE and Frozen-sample for both tumor (T) and normal (N) sample sets. Total number of present calls is on the bottom, the middle number in the circle overlap corresponds to concordant calls, and present calls unique to each subgroup are in the circles above its respective name.

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
그림 2. Venn diagram of concordance of differentially expressed genes between FFPE and Frozen samples. Genes are called differentially expressed on Affymetrix HGU133A 2.0 between tumor and normal in the FFPE and Frozen samples if they meet the following criteria: 2x or greater RMA expression value difference and p-values lower than 0.001 from a pair wise t-test. There were only 4 "false positives" in the FFPE comparisons (1 in T>N, 3 in N>T).

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
그림 3. Ingenuity pathway analysis of differentially regulated canonical pathways (A) and cellular functions (B) in tumor and normal tissues in Frozen and FFPE datasets. Pathways and cellular functions are listed along the X-axis. The Y-axis represents the $-\log(\text{significance score})$, which is defined by Ingenuity and depends on the number of genes (components) involved in the pathway/process and their relative weight in the process. The orange line represents the significance threshold. Nine of eleven canonical pathways represented in Fresh Frozen samples were detected in the FFPE samples.

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
NuGEN's Portfolio of NGS Products...

..to enable simple, fast, and scalable NGS workflows