Ovation Prokaryotic RNA-Seq System : Whole Transcriptome Profiling in Bacteria

Chris Armour, Himanshu Sethi, Victor Sementchenko, Matt Biery and Chris Raymond, NuGEN Technologies, Inc.

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
그림 1. Bacterial genomics. NGS technologies afford the opportunity to integrate phenotypes, genome sequences and global transcriptome profiling. Genome sequencing and transcriptome profiling are tightly coupled (heavy arrow) because transcriptome analysis requires the genome sequence for sequence alignments, and conversely, the interpretation of unique genomic features is guided by global transcriptome analysis.

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Ovation Prokaryotic RNA-Seq System 특징

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
그림 2. Construction of NGS transcriptome libraries. First- and second-strand cDNA Synthesis with selective priming to enrich non-rRNA transcripts from bacterial and archaeal total RNA inputs. Following end repair, the double strand cDNA is compatible with NuGEN's Encore NGS Library Systems that allows easy sample multiplexing within a single-end, transcriptome sequencing format.


그림 3. Ribosomal RNA depletion from an NGS transcriptome library. In a total RNA sample that is 95% rRNA and 5% uniquely mapping mRNA, elimination of 90% of the rRNA results in an NGS library that is 65% rRNA and 35% unique RNAs. The key feature of this library is a 7-fold enrichment of unique RNA species.


그림 4. The Ovation Prokaryotic RNA-Seq System. (A) Performance in bacteria with low GC content in their genomes. (B) Performance in bacteria with ~50% GC content in their genomes. (C) Performance in bacteria with high GC content in their genomes. Samples labeled with SP were processed with the selective priming technology in the Ovation Prokaryotic RNA-Seq System.

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
그림 5. Uniformity of transcript coverage. The horizontal axis shows a portion of the *E. coli* genome that includes the *kduD-kdul* bicistronic operon. The vertical axis shows the density of tag reads on a logarithmic scale. The sequencing results of three independent transcriptome libraries are shown; one is the random primer control and the other two are replicates of libraries produced with the Ovation Prokaryotic RNA-Seq methodology.



그림 6. Comparison of Ovation Prokaryotic RNA-Seq System selective priming (SP) approach with MICROExpress. (A) Left: Depletion of rRNA reads with MICROExpress alone, SP alone, or sequential treatment of samples with both methods. Right: Correlation between independent samples exposed to both depletion methods. (B) Left: Venn diagram of all 4401 annotated ORFs in *E. coli* showing genes detected by both methods, one method alone, or neither method. Right: Correlation of RPKM values between the two rRNA depletion methods. The inset shows 7% of genes differ by greater than 5-fold in RPKM values and less than 2% of genes differ by greater than 10-fold.

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