## Technical Report

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Neuro-2a cells transfected with 4  $\mu$ l *TransIT*-Neural Transfection Reagent and 2  $\mu$ g pEGFP (Clontech) per well of a 6-well plate and incubated in complete media. The image on the left was acquired 48 hours post-transfection on a Zeiss Axiocam. The image on the right was acquired 48 hours post-transfection on a Zeiss LSM 510 Confocal Microscope.

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DI-TNC1 astrocytes were transfected with 1  $\mu$ g *Label IT* Tracker Cy5 labeled pEGFP and 3  $\mu$ l *TransIT*-Neural Transfection Reagent per well of a 6-well plate in complete media. The blue represents the Cy5 fluorescence, and the green represents pEGFP expression. Images were acquired at 24 hours post-transfection using a Zeiss LSM 510 Confocal Microscope.

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Efficiency of knock out in various neural lines using the *TransIT*-TKO Transfection Reagent with 5-10 nM anti-firefly luciferase siRNA duplexes (Dharmacon Research, Inc). Cell lines were first transfected with plasmids expressing firefly and sea pansy luciferase using the *TransIT*-Neural Transfection Reagent. Twenty-four hours post-transfection, cell lysates were assayed for both firefly and sea pansy luciferase expression using Promega's Dual-Luciferase Reporter Assay System.

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Human astrocytes (ScienCell Research Laboratories) were seeded on poly-D-lysine coated coverslips in a 24-well plate and incubated overnight. One  $\mu$ l of *TransIT*-TKO Transfection Reagent was diluted in 50  $\mu$ l of serum-free media and 50 nM (final concentration in well) of *Label IT* siRNA Tracker Cy 3 labeled siRNA duplex was added to the diluted *TransIT*-TKO Transfection Reagent for complex formation. Subsequently, these complexes were added to 30% confluent human astrocytes in their complete growth media. Confocal image of the single astrocyte was acquired 24 hours post-transfection on a Zeiss LSM 510 Confocal Microscope.

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C6 glioma cells were seeded on poly-D-lysine coated coverslips in a 6-well plate and incubated overnight. Eight  $\mu$ l of *TransIT*-TKO Transfection Reagent were diluted in 200  $\mu$ l of serum-free media and 50 nM (final concentration in well) of *Label IT* siRNA Tracker Cy3 labeled siRNA duplex was added to the diluted *TransIT*-TKO Transfection Reagent for complex formation. Subsequently, these complexes were added to 45% confluent C6 cells in their complete growth media. Confocal image was acquired 20 hours post-transfection on a Zeiss LSM 510 Confocal Microscope.

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### *TransIT*-Neural Transfection Reagent

Code	제품명	용량
MIR 2144	<i>TransIT</i> -Neural Transfection Reagent	0.4 ml
MIR 2140	<i>TransIT</i> -Neural Transfection Reagent	1 ml
MIR 2145	<i>TransIT</i> -Neural Transfection Reagent	1 ml×5
MIR 2146	<i>TransIT</i> -Neural Transfection Reagent	1 ml×10

### *TransIT*-TKO Transfection Reagent

Code	제품명	용량
MIR 2154	<i>TransIT</i> -TKO Transfection Reagent	0.4 ml
MIR 2150	<i>TransIT</i> -TKO Transfection Reagent	1 ml
MIR 2155	<i>TransIT</i> -TKO Transfection Reagent	1 ml×5
MIR 2156	<i>TransIT</i> -TKO Transfection Reagent	1 ml×10

### *Label IT* Tracker Intracellular Nucleic Acid Localization Kit

Code	제품명	용량
MIR 7020	<i>Label IT</i> Tracker Cy3 Kit	1 Kit
MIR 7021	<i>Label IT</i> Tracker Cy5 Kit	1 Kit
MIR 7022	<i>Label IT</i> Tracker CX-Rhodamine Kit	1 Kit
MIR 7023	<i>Label IT</i> Tracker TM-Rhodamine Kit	1 Kit
MIR 7024	<i>Label IT</i> Tracker Biotin Kit	1 Kit
MIR 7025	<i>Label IT</i> Tracker Fluorescein Kit	1 Kit

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### *Label IT* siRNA Tracker Intracellular Localization Kit

Code	제품명	용량
MIR 7212	<i>Label IT</i> siRNA Tracker Cy3 Kit without Transfection Reagent	50 $\mu$ g 용
MIR 7213	<i>Label IT</i> siRNA Tracker Cy5 Kit without Transfection Reagent	50 $\mu$ g 용
MIR 7214	<i>Label IT</i> siRNA Tracker CX-Rhodamine Kit without Transfection Reagent	50 $\mu$ g 용
MIR 7215	<i>Label IT</i> siRNA Tracker TM-Rhodamine Kit without Transfection Reagent	50 $\mu$ g 용
MIR 7216	<i>Label IT</i> siRNA Tracker Fluorescein Kit without Transfection Reagent	50 $\mu$ g 용
MIR 7217	<i>Label IT</i> siRNA Tracker Biotin Kit without Transfection Reagent	50 $\mu$ g 용

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