Technical Note

FlashGel® System for DNA Recovery

By Mary Riley and Hugh White, Lonza Rockland, Inc.

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

Procedure

1. Load samples in top tier of wells.
2. Run until samples reach the second tier of wells.
3. Stop the run and add FlashGel Recovery Buffer.
4. Extract DNA from wells by pipette

Panel A – DNA Concentration Range


1) Pre-Recovery

2) Post-Recovery

3) Recovered Samples

Panel B – DNA Size Range


그림 1. FlashGel System recovery capabilities. A: Images 1 & 2 show the FlashGel Recovery Cassette before and after recovery of serial dilutions of a 1000 bp fragment. Image 3 shows analysis of 5% of each recovered sample using a 1.2% FlashGel DNA Cassette. B: 100 ng samples of fragments ranging in size from 50 bp to 4000 bp were recovered, and 3 μ l aliquots of the recovered DNA were separated on a 1.2% FlashGel DNA Cassette, along with the FlashGel QuantLadder and FlashGel DNA Marker (100 bp to 4 kb).

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
그림 2. Spin column comparison. 300 ng of 1000 bp DNA fragment were recovered using The FlashGel Recovery System (FG) or spin column type recovery products (C1 and C2). This image shows comparison of 10% of the recovered samples separated on a 1.2% FlashGel DNA Cassette, along with the FlashGel QuantLadder and FlashGel DNA Marker (100 bp to 4 kb).

Time comparison. FlashGel® System vs. Column Recovery.


그림 3. Recovery and cloning comparison. Samples of PstI/BamHI cut pBR322 were separated and 3.2 kb DNA fragments were recovered using the FlashGel Recovery System (FG) or spin column kits (C1 and C2). A: shows the comparison between 5% of each recovered DNA sample separated on a 1.2% FlashGel DNA Cassette, along with the FlashGel QuantLadder. B: PstI/BamHI cut plasmid samples from colonies transformed with recovered DNA samples; lane "V" is digest of vector with no insert. Other lanes contain the FlashGel QuantLadder and FlashGel DNA Marker (100 bp to 4 kb).

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
그림 4. Cloning of recovered PCR product. A 500 bp amplicon was amplified and recovered using the FlashGel Recovery System. 8 μ l of recovered PCR sample was ligated into the vector and transformed into *E. coli* Competent Cells.

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Panel A – Recovered PCR Products Panel B – Reamplified PCR Products


그림 5. PCR amplification of recovered DNA. Reagents from the PCR Kit were used to amplify 300 and 500 bp fragments in a multiplex PCR amplification reaction. 6 μ l of the PCR reaction was separated on the FlashGel Recovery System and the PCR product bands were recovered. A: 4 μ l of both recovered PCR fragments. The recovered fragments were used as templates for new PCR amplifications. B: 0.5 μ l aliquots of these amplification reactions. The resulting amplification products were observed for all reactions by gel analysis on a 1.2% FlashGel DNA Cassette, along with the FlashGel QuantLadder. Reactions 1 & 3 used primers for 300 and 500 bp fragments while reactions 2 & 4 used only single primer set.