탁월한 Adenovirus 정제 Adenovirus Purification Kit

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Achieve High Titer of High Quality Virus

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
그림 2. Time course of adenoviral yields from cell pellets and supernatants, and the concomitant progression of the cytopathic effect. HEK-293 cells were seeded on a 6-well plate using 1x10⁵ cells/cm² and 0.28 ml/cm² medium, and incubated for 16 hr. They were then transduced with an Adeno-X Expression System 2 construct expressing DsRed-Express at a multiplicity of infection (MOI) of ~1-2, and incubated for 2-4 days. Panel A, The cytopathic effect was followed by examining the infected HEK 293 cells by phase contrast microscopy (left column) and expression was monitored by DsRed-Express fluorescence (right column) on days 2-4. It is clear that by day 4, the cytopathic effect was nearly complete. Panel B, On each of days 2-4, the cells were pelleted, and their pellets and supernatants collected. The pellets were resuspended in fresh medium, lysed by freeze-thawing three times to release the viral particles, and then centrifuged. The viral titer (total IFU) from each pellet and supernatant was determined by infecting fresh HEK-293 cells, and counting the cells under a fluorescence microscope 48 hr following infection. The adenoviral yields from each of the pellets and supernatants from days 2-4 were then compared. Most of the virus was found in the cell pellet regardless of when it was collected. Maximal amounts of virus were found in the cell pellet from day 4, at which time the supernatant contained only ~2% of that found in the pellet.

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
그림 3. Mass balance analysis of our new virus purification process. HEK-293 cells were seeded and infected with Adeno-XTM DsRed-Express as described in 그림 2, and harvested when the cytopathic effect was nearly complete. The resulting adenovirus was then purified according to the Adeno-XTM Maxi Purification Kit protocol. Total IFU was determined (as described in 그림 2) for each fraction (crude, flowthrough, wash, and eluate), and plotted. 72% of the starting material was recovered in the eluate as purified virus. The process was efficient and fast, allowing us to purify 4.68 x 10¹⁰ adenoviral particles in one hour.

Product	Size	TaKaRa Code
Adeno-X TM Maxi Purification Kit	2 preps	631532
	6 preps	631533
Adeno-X TM Expression System 2	5 rxns	631524
Adeno-X TM ViraTrak DsRed-Express Expression System 2	10 rxns	632516

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