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
그림 1. Ligand-dependent, targeted, and reversible protein stabilization. A small destabilization domain (DD) is fused to a target protein of interest. The small membrane-permeable ligand Shield1 binds to the DD and protects it from proteasomal degradation. Removal of Shield1, however, causes rapid degradation of the entire fusion protein. The default pathway for the ProteoTuner Systems is degradation of the fusion protein, unless Shield1 is present.

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Protein Regulation	Technology	CLONTECH Products
Transcriptional	Controllable Gene Expression	Tet-On®/Off® Systems 36-48 hr
Translational	RNA Interference	Knockout™ (Inducible) RNAi Systems 48-72 hr
Post-translational	Ligand-Dependent Protein Stabilization / Degradation	ProteoTuner™ 30 min-1 hr

Protein Regulation Technologies Available from CLONTECH

그림 2. ProteoTuner™ System is more rapid and targeted control of protein expression compared to transcriptional or translational regulation.

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
그림 3. The fluorescence intensity of DD-DsRed-Express is directly related to the concentration of the stabilizing ligand Shield1. Cells were infected with pRetroX-PTuner DsRed-Express and treated with different concentrations of Shield1. The amount of DD-tagged DsRed-Express stabilized by different concentrations of Shield1 was detected by fluorescence intensity 18 hr later, using a BD FACSCalibur™ flow cytometer. RFU = relative fluorescence units.


그림 4, DD-DsRed-Express fluorescence is directly related to the concentration of the stabilizing ligand Shield1. In order to visualize the amount of DD-DsRed-Express expressed in a cell, cells were infected with pRetroX-PTuner DsRed-Express IRES ZsGreen1 and treated with different concentrations of Shield1. The amount of DD-tagged DsRed-Express stabilized by different concentrations of Shield1 was detected via Western blot using the Living Colors DsRed Polyclonal Antibody. Lane 1 : molecular weight marker, Lane 2 : 1X loading buffer, Lane 3: untreated HeLa cells (no virus, no Shield1). Lane 4 : HeLa cells infected with the DD-DsRed Express construct; no Shield1, Lanes 5,8 : HeLa cells infected with the DD-DsRed-Express construct and treated with 50, 250, 500, and 1,000 nM Shield1 respectively. Lane 9 : 1X loading buffer, Lane 10 : HEK 293 DsRed-Express stable cell line.

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¹ Department of Chemistry, Stanford University, Stanford, California 94305, USA

² Department of Chemical and Systems Biology, Stanford University, Stanford, California 94305, USA

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
그림 5. FKBP mutants display Shield1-induced stability. Cells expressing FKBP-mutant YFP fusions and a negative control were either mock-treated or treated with 1 μM Shield1 for 24 hr. Lysate was analyzed by Western blot with an anti-FKBP antibody. Differing degrees of Shield1-induced stabilization could be observed for the different mutants analyzed. FKBP-YFP was not detected in lysate from mock-treated cells.


그림 6. Kinetics of Shield1-dependent stabilization. Panel A, NIH 3T3 cells stably expressing FKBP-mutant YFP fusions were treated with threefold dilutions of Shield1 (1 μM to 0.1 nM). Mutant V24A showed the most efficient rescue (EC₅₀ ~5 nM), whereas the more destabilizing L106P required higher concentrations of Shield1 (EC₅₀ ~100 nM) to stabilize the YFP fusion protein. Panel B, NIH 3T3 cells stably expressing FKBP-YFP fusions were treated with 1 μM Shield1 for 24 hr and then washed with media to remove Shield1. L106P was degraded most quickly, with protein levels becoming negligible within 4 hours. Changes in fluorescence were monitored by flow cytometry. MFI = mean fluorescence intensity.

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
그림 7. Comparison of RNAi and the ProteoTuner systems. The time required for RNAi knockdown of lamin A/C (Panel A) was compared to the time required for degradation of DD-YFP upon removal of Shield1 from NIH 3T3 cells stably expressing the fusion protein (Panel B).

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
그림 8. Predictable and reversible, Shield1-dependent regulation of intracellular protein levels. NIH 3T3 cells stably expressing DD-YFP were treated with varying concentrations of Shield1 over the course of one week, and samples of the population were assayed by flow cytometry at the indicated time points. Predicted fluorescence is based upon the dose response experiment shown in Figure 2A. MFI = mean fluorescence intensity.

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
그림 9. Stabilization of DD-tagged Cdc42 (Q61L) with Shield1 results in predictable changes in cellular morphology. NIH 3T3 cells stably expressing fusions of a constitutively active small GTPase to the DD were split into three pools. The first population (Panel A) was mock-treated and the second population (Panel B) was treated with 1 μM Shield1 for 24 hr. The third population (Panel C) was treated with 1 μM Shield1 for 24 hr, then washed with media and cultured in the absence of Shield1 for 48 hr. Cells were serum-starved for 12 hr, fixed, stained with Alexa Fluor® 488-conjugated phalloidin, and visualized using confocal microscopy.

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
그림 10. Functional analysis of essential proteins using the ProteoTuner System. Immunofluorescence and Western blot analyses were performed on *T. gondii* parasite strains expressing DD-MyoA (Panel A) and DD-Rab11ADN (Panel B). Both visualization methods confirm that the DD fusion proteins are stabilized in the presence of Shield1. Immunofluorescence analyses were performed 4 hr after induction with Shield1, using monoclonal anti-myc (clone 9E10, Sigma) and either anti-Gap45⁹ or polyclonal anti-IMC1¹⁰. Scale bar, 10 μm. Western blot analyses using these antibodies were performed on the same parasite strains, grown in the presence or absence of Shield1 for 48 hr (DD-MyoA) or for two hr (DD-Rab11ADN) before analysis.

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ProteoTuner System - Technical Note 2

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²Department of Chemistry and Systems Biology, Stanford University, 269 Campus Drive, Stanford, California 94305, USA.

³These authors contributed equally to this work.


그림 11. Quantification of invasion and replication of parasite strain DD-Rab11ADN. Total number of parasitophorous vacuoles was determined. Mean values of three independent experiments ± s.d. are shown (?/+, parasites treated with Shield1 after invasion; +/-, parasites treated with Shield1 before and after invasion). Asterisks indicate a significant difference in total invasion between parasites not treated with Shield1 and parasites treated with Shield1 before and after invasion (P < 0.01, two tailed Student's t-test).

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표 1 : Six Single-Vector ProteoTuner™ Systems Now Available

System	TaKaRa Code	Vector Type	Antibiotic Resistance ¹	Fluorescent Protein
ProteoTuner™ System	632172	Plasmid	Kanamycin /G418	None
ProteoTuner™ IRES2 System	632168	Plasmid	Kanamycin /G418	AcGFP1
Retro-X™ ProteoTuner™ System	632171	Retroviral	Ampicillin /Puromycin	None
Retro-X™ ProteoTuner™ IRES System	632167	Retroviral	Ampicillin /None	ZsGreen1
Lenti-X™ ProteoTuner™ System	632173	Lentiviral	Ampicillin /Puromycin	None
Lenti-X™ ProteoTuner™ Green System	632175	Lentiviral	Ampicillin /None	ZsGreen1

¹ Bacterial/eukaryotic

ProteoTuner™ System Components (TaKaRa Code 632172)

- pPTuner Vector
- Shield1

ProteoTuner™ IRES2 System Components (TaKaRa Code 632168)

- pPTuner IRES2 Vector
- Shield1

Retro-X™ ProteoTuner™ System Components (TaKaRa Code 632171)

- pRetroX-PTuner Vector
- Shield1

Retro-X™ ProteoTuner™ IRES System Components (TaKaRa Code 632167)

- pRetroX-PTuner IRES Vector
- Shield1

Lenti-X™ ProteoTuner™ System (TaKaRa Code 632173)

- pLVX-PTuner Vector
- Shield1

Lenti-X™ ProteoTuner™ Green System (TaKaRa Code 632175)

- pLVX-PTuner Green
- Shield1

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제품명	Size	TaKaRa Code
Shield1	60 µl	631037
	200 µl	631038
Living Colors® DsRed Polyclonal Antibody	100 µl	632496
Tet-On® Advanced Inducible Gene Expression System	each	630930
Tet-Off® Advanced Inducible Gene Expression System	each	630934
Knockout™ Single Vector Inducible RNAi System	each	630933

※ License Notice [5] [7] [8]

Real Time Gene Expression Measurement in Living Cells

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ATTO AB-2550 Kronos Dio