continued...

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표 2: Surface Protein Phenotypes of Some Human Stem Cells

Gene	hCD34 ⁺	hMSC	hADSC
VLA-4 ¹	+	-	+
VLA-5 ¹	+	+	+
Pit-1 ²	0.36	1	0.31

¹ VLA-4 and VLA-5 are integrins that bind RetroNectin.

² Pit-1 is the cellular receptor for the GaLV Env protein. These values represent Pit-1 expression relative to that of hMSC.


그림 2. RetroNectin-mediated transduction outperforms standard methods of transduction. A GaLV-pseudotyped retrovirus carrying a fluorescent protein marker was used to transduce human hematopoietic (hCD34⁺), mesenchymal (hMSC), and adipose (hADSC) stem cells, all of which express different levels of Pit-1 (the GaLV receptor) (see Table II). In each case, the RBV method transduced the stem cells more efficiently than either the polybrene- or the protamine-based method.

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제품명	용량	TaKaRa Code
RetroNectin Reagent	0.5 mg	T100A
	2.5 mg	T100B
RetroNectin Precoated Dishes	10 dishes	T110A

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- Hanenberg, H. et al. (1996) *Nature Med.* 2(8):876 – 882.
- Hanenberg, H. et al. (1997) *Hum. Gene Ther.* 8(1):2193 – 2206.
- Chono, H. et al. (2001) *J. Biochem.* 130(3):331 – 334.

* License Notice: [16]


그림 3. Effective retroviral and lentiviral transduction of Jurkat cells on RetroNectin-coated plates. Jurkat cells were transduced with either retrovirus (Panel A) or lentivirus (Panel B) by using serial RBV-based transductions. Each VSV-G pseudotyped virus was engineered to express ZsGreen1. Viruses were bound to RetroNectin-coated plates for 4 hr followed by incubation with the cells for 24 hr. Following the first infection, the cells were collected and added to a second RetroNectin plate containing freshly bound virus. At 96 or 72 hr after the second round of retroviral or lentiviral infection, respectively, the cells were analyzed for fluorescent protein expression using flow cytometry. The transduced cells shown in Panel A were 75.5% positive with an MFI of 138; those shown in Panel B were 94.8% positive with an MFI of 120.