# Lentiviral delivery system을 위한 total solution

## Lenti-X™ Expression System

### ■ 고효율의 Lentiviral Packaging System

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
그림 1. The Lenti-X HT Packaging System. A lentiviral vector (e.g. pLVX-Puro) and the Lenti-X HT Packaging Mix are cotransfected into 293T cells using the highly efficient LentiPhos HT transfection system. High titer supernatants are ready for use 48 hr after transfection.

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
그림 2. Lenti-X vectors for many applications. Lenti-X vectors contain sequence elements that facilitate lentiviral packaging, boost transgene expression, or both. Among them are the HIV-1 LTRs and packaging signal (Ψ), a Rev response element (RRE), the central poly-purine tract/central termination sequence (cPPT/CTS), and the woodchuck hepatitis virus posttranscriptional regulatory element (WPRE). All vectors are designed to be used with the Lenti-X HT Packaging System and the Lenti-X 293T Cell Line, which will produce high titers of VSV-G-pseudotyped lentivirus for transducing virtually any cell type. See article text for vector descriptions and applications.

# Lentiviral delivery system을 위한 total solution Lenti-X™ Expression System

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
그림 3. High infectivity of supernatants produced by the Lenti-X HT Packaging System. The Lenti-X HT Packaging System (Panel A) and a packaging system from a competitor (Panel B) were each used to generate viral supernatants from their respective lentiviral system vector that was engineered to express the ZsGreen1 fluorescent protein. As little as 10 μl of supernatant from the Lenti-X HT Packaging System transduced the majority of these HeLa cells, whereas 10 μl of supernatant from the other system transduced only a small percentage of the cells. Transduced cells were quantified by flow cytometry


그림 4. Transduction of neural progenitor cells by Lenti-X Lentivirus. Recombinant lentivirus for expressing ZsGreen1 was produced from the Lenti-X HT Packaging System and used to transduce normal human neural progenitor cells. A single transduced cell is shown under phase contrast microscopy (Panel A) and fluorescence microscopy (Panel B).

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
그림 5. Induced expression in the Tet-Off Advanced and Tet-On Advanced Systems. The Tet-controlled transactivators are fusion proteins that contain a DNA-binding TetR domain joined to three minimal transcription activation domains from HSV VP16. Each transactivator has been optimized for expression in mammalian cells. In Tet-Off Advanced Systems, the basal state is maintained in the presence of doxycycline (Dox), and induced by its withdrawal. Tet-On Advanced Systems are activated in the presence of Dox. System induction produces high-level transcription of your gene from PTight.

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### ■ Fluorescent Protein Coexpression & Tag

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

Figure 1. The Lenti-X GoStix protocol takes only 30 seconds.

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