Lentivirus Systems & Tools

Lentivirus Systems

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그립 1. 4th generation lentiviral packaging system. A lentiviral vector and the Lenti-X HTX Packaging Mix are cotransfected into 293T cells. High titer lentiviral supernatants are ready for use 48 hr after transfection.

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그립 2. High infectivity of supernatants produced by Lenti-X. Lenti-X (Panel A) and a packaging system from a competitor (Panel B) were each used to generate virus containing a vector system for expressing the ZsGreen1 fluorescent protein. 10 µl of supernatant from each system was used to transduce HeLa cells. ZsGreen1-positive cells were quantified by flow cytometry Lenti-X transduced the majority of cells, whereas the other system transduced only a small percentage of the cells.

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그림 4. Efficient concentration with minimal loss. Lentiviral supernatant from a pLVX–ZsGreen1 vector was concentrated from 3 ml down to 30 µl using the Lenti–X Concentrator reagent, which reflected a 100–fold increase in viral titer. Measuring the total amount of virus contained in each sample indicated that the resuspended pellet captured 90% of the virus present in the original sample. Samples were titrated using HT1080 cells and analyzed by flow cytometry 48 hr post–transduction.

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표1. Lenti-X Concentrator vs. Ultracentrifugation

Feature	Lenti-X Concentrator	Ultracentrifugation
Easily Scalable	Yes	No
Specialized Equipment	No	Yes
Time Required	~1 hr	4 hr to overnight
Ease-of-Use	++++	+
Yield	>90%	>90%

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그림 5. The Lenti-X Maxi Purification Kit allows you to generate high yields of purified lentivirus from crude packaging cell supernatants. The gravity column-based method (bind, wash, elute) is extremely simpleand effective, and preserves virus infectivity much better than filter-based methods.

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그림 6. Lenti-X Maxi Purification Kit yields are far higherthan the yields of filter-based methods. Panel A. Virus contentin the indicated Lenti-X column fractions was tracked using either Lenti-X qRT-PCR (RNA copies) or flow cytometry/fluorescence (IFU) titration. The mean values from seven experiments are shown. Panel B. In a headto-head comparison, Lenti-X column purification recovers more virus than a filter-based method.far fewer contaminating proteins at equivalent IFU than either sample prepared from the filter-based systems

그립 7. The Lenti-X Maxi Purification Kit yields highly purified lentivirus. Equivalent amounts of purified virus (1 x 10⁵ IFU) prepared using either the Lenti-X Maxi Purification Kit (Lenti-X) or a filter-based system (Vendor M & Vendor S) were subjected to SDS-PAGE and silver-stained. The Lenti-X sample clearly contained far fewer contaminating proteins at equivalent IFU than either sample prepared from the filter-based systems.

Instantly Test for Lentivirus Lenti–X™ GoStix™

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그림 8. p24 is located in the lentiviral capsid and is one of 4 proteins encoded by the HIV-1 gag gene.

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그립 9. Higher proviral content is associated with higher expression levels. Provirus content in the genomic DNAs of HT1080 cells transduced at either high or low MOI were determined using the Lenti–X Provirus Quantitation Kit. AcGFP1 expression levels (MFI) were determined using flow cytometry. The figure represents data pooled from two groups having either low-copynumber ($\langle 3; n = 8 \rangle$ or high-copy-number ($\rangle 9; n = 7$) proviruses.

그립 10. Equivalent MOI in different cell types can result in different provirus copy numbers. HT1080, HeLa, and HEK 293 cell cultures were transduced at equivalent MOIs using a single lentiviral stock (100 µl). The provirus copy numbers in the genomic DNAs of the transduced cell populations were determined using the Lenti-X[™] Provirus Quantitation Kit.

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그립 11. Proviruses amplified by inverted PCR correspond to qPCR copy number. Provirus copy numbers were determined for genomic DNA purified from 106 stably transduced HT1080 cells representing either individual clones (Lanes 1-3) or polyclonal populations (Lanes 4 & 5). The DNA samples were also subjected to inverted PCR analysis to amplify individual provirus junction sequences using vector-specific primers (1). NTC = no template control.

Note: Not all proviruses can be amplified with this method.

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Magnetic Separator

그림 12. Lenti-X Accelerator provides high transduction efficiency in a 25 min protocol. Lentiviral transduction of HT1080 cells was carried out for 5 min with Lenti-X™ Accelerator beads after a 20 min incubation to bind the beads to the virus—and for 5 min or overnight with Polybrene. After the cultures were grown for an additional 72 hr at 37°C, the number of transduced cells was determined by flow cytometry.

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Viral Receptor Boosters are concentrated exosome-like vesicles that are applied to target cells prior to infection with virus. Booster treatment increases the cell surface density of the receptor recognized by the infecting virus, thus increasing transduction efficiency. Using this technology, for example, you can coat human cells with the ecotropic receptor (which is otherwise absent), enabling them to be transduced with ecotropic retrovirus or ecotropic lentivirus.

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그립 15. Viral transduction of human cells with ecotropic lentivirus following Ecotropic Receptor Booster treatment. HT1080 cells were seeded in 6-well plates 24 hr prior to transduction and incubated with 10 µl Ecotropic Receptor Booster for 2 hr. Cells were then transduced with Lenti-X ZsGreen1 lentivirus (MOI=~15) produced using Clontech's Lenti-X HTX Ecotropic Packaging System (Code 631251) and assayed 48 hr later for ZsGreen1 expression.

그립 16. Ecotropic Receptor Booster aides transduction of ecotropic lentivirus. Panel A. Usually, it is not possible to infect Jurkat cells, hMSCs, or NHNPs with ecotropic pseudotyped lentivirus. However, pre-treatment of these cell types with Ecotropic Receptor Booster allows for very efficient transduction. Panel B. A panel of human cell lines, normally resistant to transduction by ecotropicpseudotyped lentiviral vectors, were efficiently transduced after treatment with Ecotropic Receptor Booster.

Lenti-X[™] Integration Site Analysis Kit

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