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
그림 2. Constructing recombinant adenovirus with In-Fusion technology. DNA sequences can be rapidly transferred as PCR products to any pAdenoX vector using the In-Fusion cloning method. In this example, your gene of interest is amplified with 15 bp extensions that are homologous to the ends of the linearized adenoviral vector. The PCR product is then purified and mixed with the linearized adenoviral vector of choice in the In-Fusion reaction. Following the reaction, a portion of the mixture is transformed into Stellar Competent Cells and screened. Once a PCR-positive clone is identified, the recombinant pAdenoX vector is amplified, purified, and subsequently linearized with the restriction enzyme PaclI, then transfected into HEK 293 cells for viral rescue and amplification.

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
표 1 : Adeno-X Adenoviral System 3 Formats

Code	Product	Description	Vector Map
631180	Adeno-X Adenoviral System 3 (Tet-On 3G Inducible)	<ul style="list-style-type: none"> Tightly-controlled, doxycycline-inducible expression system 	<p>Tet-On 3G Transactivator</p>
632269	Adeno-X Adenoviral System 3 (CMV)	<ul style="list-style-type: none"> Constitutive expression from a CMV promoter 	
632268	Adeno-X Adenoviral System 3 (CMV, Red)	<ul style="list-style-type: none"> Constitutive expression from a CMV promoter Red fluorescent protein to easily monitor transfection and transduction 	<p>DsRed-Express</p>
632267	Adeno-X Adenoviral System 3 (CMV, Green)	<ul style="list-style-type: none"> Constitutive expression from a CMV promoter Green fluorescent protein to easily monitor transfection and transduction 	<p>ZsGreen1</p>
632266	Adeno-X Adenoviral System 3 (Universal)	<ul style="list-style-type: none"> Use any promoter and any polyA sequence Ideal for tissue-specific expression or expression of shRNA or miRNA 	
632265	Adeno-X Adenoviral System 3 (Universal, Red)	<ul style="list-style-type: none"> Use any promoter and any polyA sequence Ideal for tissue-specific expression or expression of shRNA or miRNA Red fluorescent protein to easily monitor transfection and transduction 	<p>DsRed-Express</p>
632264	Adeno-X Adenoviral System 3 (Universal, Green)	<ul style="list-style-type: none"> Use any promoter and any polyA sequence Ideal for tissue-specific expression or expression of shRNA or miRNA Green fluorescent protein to easily monitor transfection and transduction 	<p>ZsGreen1</p>

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그림 3. The Adeno-X Tet-On 3G Systems generate very highfold induction, with up to 3000-fold difference between induced and uninduced states. Using equal amounts of high-titer supernatants, HeLa cells cultured at the indicated concentrations of Dox were infected with Adeno-X Tet-On 3G Luciferase virus. Cultures were harvested and assayed for luciferase activity.

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
그림 4. The Universal Adeno-X Expression Systems contain vectors that lack a promoter and polyA signal in the cloning site. You can either clone an expression cassette from a preexisting construct into the vector (Panel A) or create a new one using multiple fragment In-Fusion HD cloning (Panel B).

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표 2 : Multiple Fragment Cloning Efficiency with the Adeno-X System 3

	Number of colonies (1/10th plated)	% correct
pAdenoX + 1 insert	200-300	90
pAdenoX + 2 inserts	200-300	60
pAdenoX + 3 inserts	25-40	25

표 3 : Adenoviral Gene Delivery vs. Lentiviral Gene Delivery

	Lentivirus	Adenovirus
Infects many different human cell types	Yes	Yes
Infects both dividing and non-dividing cells	Yes	Yes
Non-integrating virus	No	Yes
High level of protein expression (up to 10-20% total protein)	No	Yes
Ability to accommodate long inserts (up to 8 kb)	No	Yes
Easy to scale-up/amplify	No	Yes
Easy to get titers >10 ⁹ IFU/ml	No	Yes
Easy to get a multiplicity of infection >25 copies per cell	No	Yes

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Code	제품명	용량
631180	Adeno-X Adenoviral System 3 (Tet-On 3G Inducible)	Set
632269	Adeno-X Adenoviral System 3 (CMV)	Set
632268	Adeno-X Adenoviral System 3 (CMV, Red)	Set
632267	Adeno-X Adenoviral System 3 (CMV, Green)	Set
632266	Adeno-X Adenoviral System 3 (Universal)	Set
632265	Adeno-X Adenoviral System 3 (Universal, Red)	Set
632264	Adeno-X Adenoviral System 3 (Universal, Green)	Set

License Notice [K2, K20, K45]